



DYOMICS
Colours for Life

Catalogue

Fluorescent Dyes for Bioanalytical
and Hightech Applications

8th Edition – spring 2017

Content

Mission statement	2
Labels	
– Table with spectral characteristics of DY-Labels	4
– UV-excitation	8
– Blue-green excitation	10
– Green excitation	16
– Orange excitation	24
– Red excitation	28
– Far red excitation	38
– NIR-excitation	44
– MegaStokes-Dyes for multi-colour-detection	62
– Quencher	72
PEGylated Fluorescent Labels	78
Special Dyes	
– Special Dyes	86
– Direct from Research	89
– Sensor Dyes	91
– Laser Dyes	92
Haptens	
– Fluorescently labeled biotins	93
– Fluorescently labeled phalloidins	93
Fluorescently labeled proteins and nucleotides	
– Fluorescently labeled streptavidin	94
– Fluorescently labeled antibodies for secondary detection	94
– Others	94
– Fluorescently labeled nucleotides	94
Protocols	95
Services	98
Recommended literature	99

Preface to the 8th Edition

For nearly 20 years Dyomics has successfully advanced and brought to market decades of academic learning and research in the field of organic dyes.

Decades ago xanthylium-dyes (rhodamines) and indocyanines (CY3, CY5 and CY7) were the preferred chromophores used in Life Sciences. Each of these chromophores had clear strengths, but they also had significant limitations: xanthylium-derivatives exhibit high quantum yields and photostability, and indocyanines can reach wavelengths up to the NIR region very easily. However, xanthylium-dyes suffer from a low efficiency for excitation and indocyanines show poor stability and restricted opportunities for molecular engineering. Through the design, synthesis and tailoring of innovative fluorophores, Dyomics brought together the advantages of both xanthylium and indocyanines in a new family of hybrids between xanthylium dyes and indocyanines, called benzopyrylium based fluorescent labels.



A second notable innovation was Dyomics development of the valuable MegaStokes™ dyes, which allow multiplex applications.

According to the needs of our partners and their platforms, it is possible to tailor the new DY-dyes to fulfill given demands in terms of water solubility, reduce nonspecific mutual reactions, increase thermic and photochemical stability and provide opportunities to react selectively with the desired biomolecules.

These new stable DY-dyes work up to the NIR and red region and have gained broad acceptance from customers worldwide.

Established applications can be found in cutting edge biological research and medical diagnostics including DNA-sequencing, high resolution microscopy, flow-cytometry, DNA- and/or protein labelling and analysis, and in vivo Imaging.

Dyomics rich heritage of applied research and break-through innovations continues to this day, and we passionately serve both the academic and industrial communities with our tailored DY-dyes.

Dyomics is an ISO 9001 certified, R&D focused, customer-centered company and we are look forward to providing you with DY-dyes that exceed your most exacting requirements.

Sincerely yours,

A handwritten signature in dark ink, appearing to read 'Peter Czerney'.

Dr. Peter Czerney
(Managing Director)

Mission Statement

Development and growth in the use of organic dyes in biotechnology, biomedical research and clinical diagnostics is continuing. Powerful new markers and/or probes which have been constructed using knowledge of both organic chemistry and molecular biology, are combined with automated imaging workstations to define the content, activity, and dynamics of different materials in life science. These innovative applications have led to a new dye chemistry. At the present the design of functional chromophores and fluorophores is the most significant and enticing field in dye chemistry.

With more than three decades of experience in synthesis of functional fluorescent dyes, the staff of Dyomics has an excellent scientific foundation in the search for labels for bioanalytics and diagnostics. Dyomics accepts the challenge to develop new high quality "tailor made" dyes for bioanalysis and related fields of high technology. This is, however, not the only goal of our company.

As a established company, we are open-minded and flexible in nearly every aspect of colour chemistry. Our product portfolio ranges from different labels and stains for optical sensing, which includes highly hydrophilic probes and dyes including new MegaStokes' dyes with enhanced Stokes' shift.

In cooperation with business partners in the area of molecular biology and pharmaceutical research we continue the development of new tailored products. Customers' requests for special modifications and/or synthesis are always welcome and belong to our day by day business.

Our greatest wish is that our products and services will satisfy you as the customer, so that we may establish a lasting relationship with each one of you.

Products

The production of new chromophores has extensively changed the possibilities for detection of biomolecules and their interaction in recent years and has enabled previously impossible insights to be made.

The advantages of optical detection methods are in the speed of the measuring process and the high sensitivity (up to single molecule detection). In principal, the collection of data over time as well as parallel or multiple readings are feasible. Hence it is possible to answer an array of questions in biotechnology or medical diagnostics in a cost-saving and reliable manner including the application of miniaturized analysers.

In this area Dyomics offers technical know how for diverse users and hardware developers in the choice and design of suitable chromophores optimized for particular excitation and emission wavelengths.

Currently, there are more than 60 fluorophores and quenchers of the DY family available as NHS- and TFP-esters, maleimides and others for the labeling of biomolecules.

In particular, our company has great experience in the development of fluorophores which can be excited with low-cost diode lasers and show emission in the red and near infrared areas of the spectrum, as well as with custom made FRET systems.

According to the demands from the customer, the activated DY dyes can be coupled covalently to amino acids, antibodies, hormones, sugars, dNTPs, oligonucleotides or other affinity tags to corresponding probes or markers.

The proprietary DY dyes are especially useful in the analysis of enzyme kinetics and DNA sequences, in the evaluation of chromosomes by Fluorescence-In-Situ-Hybridisation (FISH), in Fluorescence Correlation Spectroscopy (FCS), in FRET assays, in microscopy and in DNA and protein chips although new applications are always being discovered. Current developments for example include cell organelle specific stains and probes.

Quality Statement

Dyomics is committed to a high standard of quality. Since the second year after the foundation of the company in 1999, Dyomics has worked in accordance to the quality management system DIN EN 9001:2000 (www.tuevclub.com)

The continuous monitoring and documentation of all processes, the exclusive use of certified reagents and the final control of the products using ^1H - and ^{13}C -NMR, LC/MS, HPLC, UV-VIS and fluorescence spectroscopy guarantee a high conformity and quality of every single product, a guarantee for reproducible results.

Quality management relates not only to the different products, but also to the rapid processing of both routine and R&D contracts. Characteristic therefore is the obliging and customer orientated service and the efficient implementation of synthesis. Ultimately the key to the fulfilment of our quality claims is our qualified and experienced work force.



Table with spectral characteristics of DY-Labels

label	absorption maximum* [nm]	molar absorbance* [l·mol ⁻¹ ·cm ⁻¹]	emission maximum* [nm]	solvent	CF		page
					[280nm]	[260nm]	
DY-350	353	19,000	432	PBS	0,085	0,257	9
DY-351	338	33,000	437	PBS	0,302	0,405	9
DY-405	400	32,000	423	PBS	0,772	0,252	11
DY-410	405	34,000	460	PBS	0,227	0,270	11
DY-415	418	40,000	463	EtOH	0,217	0,315	12
DY-430	444	37,000	487	PBS	0,278	0,456	12
DY-431	442	35,000	496	PBS	0,247	0,479	13
DY-478	478	90,000	518	EtOH	0,186	0,216	13
DY-490	491	73,000	515	PBS	0,113	0,198	14
DY-495	494	78,000	521	pH9 buffer	0,178	0,261	14
DY-505	507	80,000	528	EtOH	0,165	0,236	15
DY-530	533	100,000	554	PBS	0,148	0,295	17
DY-546	549	100,000	569	PBS	0,213	0,266	17
DY-547P1	551	150,000	565	PBS	0,075	0,029	18
DY-548P1	554	150,000	568	PBS	0,074	0,028	18
DY-549P1	556	150,000	570	PBS	0,070	0,035	19
DY-550	558	120,000	577	EtOH	0,044	0,147	19
DY-554	544	100,000	570	EtOH	0,233	0,383	20
DY-555	547	100,000	573	EtOH	0,185	0,341	20
DY-556	548	100,000	574	EtOH	0,164	0,318	21
DY-557	556	100,000	576	EtOH	0,175	0,331	21
DY-560	560	120,000	578	EtOH	0,189	0,337	22
DY-580	572	100,000	589	EtOH	0,151	0,245	22
DY-585	584	120,000	604	EtOH	0,192	0,509	23
DY-590	581	120,000	600	EtOH	0,170	0,224	25
DY-591	581	120,000	598	EtOH	0,143	0,189	25
DY-594	594	92,000	615	PBS	0,622	0,347	26
DY-605	600	110,000	624	PBS	0,552	0,237	26
DY-610	610	80,000	632	EtOH	0,210	0,200	27
DY-615	623	200,000	643	EtOH	0,058	0,070	27
DY-630	638	200,000	658	EtOH	0,092	0,092	29
DY-631	637	200,000	657	EtOH	0,312	0,259	29
DY-632	636	200,000	658	EtOH	0,108	0,097	30
DY-633	638	200,000	658	EtOH	0,088	0,086	30
DY-634	636	200,000	657	EtOH	0,095	0,075	31
DY-635	648	200,000	670	EtOH	0,181	0,154	31
DY-636	647	200,000	670	EtOH	0,084	0,045	32
DY-641	624	80,000	654	EtOH	0,262	0,221	32
DY-643	626	80,000	646	EtOH	0,238	0,227	33
DY-647P1	648	250,000	666	PBS	0,024	0,019	33
DY-648P1	651	250,000	669	PBS	0,024	0,021	34
DY-649P1	654	250,000	672	PBS	0,025	0,018	34
DY-650	656	220,000	676	EtOH	0,058	0,113	35
DY-651	655	220,000	677	EtOH	0,052	0,101	35
DY-652	653	220,000	676	EtOH	0,060	0,110	36
DY-654	653	220,000	677	EtOH	0,067	0,145	36
DY-660P1	663	200,000	682	PBS	0,070	0,071	37
DY-675	675	180,000	699	EtOH	0,080	0,151	39
DY-676	675	180,000	699	EtOH	0,078	0,147	39
DY-677	674	180,000	698	EtOH	0,087	0,166	40
DY-678	674	180,000	694	EtOH	0,100	0,179	40
DY-679P1	679	200,000	697	PBS	0,115	0,107	41
DY-680	691	140,000	709	EtOH	0,099	0,142	41
DY-681	692	140,000	709	EtOH	0,089	0,125	42
DY-682	692	140,000	709	EtOH	0,085	0,124	42
DY-684	685	140,000	699	EtOH	0,121	0,147	43

Table with spectral characteristics of DY-Labels

label	absorption maximum* [nm]	molar absorbance* [l · mol ⁻¹ · cm ⁻¹]	emission maximum* [nm]	solvent	CF		page
					[280nm]	[260nm]	
DY-700	707	140,000	728	EtOH	0,103	0,172	45
DY-701	709	140,000	730	EtOH	0,127	0,197	45
DY-702	709	140,000	728	EtOH	0,076	0,128	46
DY-703	705	140,000	721	EtOH	0,096	0,140	46
DY-704	706	140,000	721	EtOH	0,109	0,135	47
DY-705	710	140,000	732	EtOH	0,072	0,123	47
DY-706	710	140,000	733	EtOH	0,077	0,124	48
DY-720	726	140,000	744	EtOH	0,095	0,160	48
DY-730	734	240,000	755	EtOH	0,047	0,054	49
DY-731	736	240,000	755	EtOH	0,052	0,057	49
DY-732	735	240,000	756	EtOH	0,064	0,071	50
DY-734	733	240,000	755	EtOH	0,045	0,051	50
DY-736	743	240,000	764	EtOH	0,071	0,051	51
DY-747P1	747	240,000	769	PBS	0,037	0,026	51
DY-749P1	747	240,000	769	PBS	0,037	0,026	52
DY-750	751	270,000	774	EtOH	0,030	0,087	52
DY-751	752	270,000	772	EtOH	0,029	0,088	53
DY-752	750	270,000	771	EtOH	0,028	0,080	53
DY-754	748	270,000	771	EtOH	0,035	0,106	54
DY-765	765	130,000	785	EtOH	0,218	0,304	54
DY-776	772	240,000	787	EtOH	0,074	0,144	55
DY-777	770	240,000	788	EtOH	0,070	0,138	55
DY-778	767	240,000	787	EtOH	0,081	0,149	56
DY-780	783	170,000	799	EtOH	0,079	0,082	56
DY-781	784	170,000	796	EtOH	0,075	0,074	57
DY-782	785	170,000	794	EtOH	0,096	0,090	57
DY-784	779	170,000	780	EtOH	0,102	0,110	58
DY-800	777	280,000	791	EtOH	0,045	0,034	58
DY-805	800	150,000	824	EtOH	0,141	0,107	59
DY-820	823	200,000	843	EtOH	0,090	0,064	59
DY-831	844	220,000	875	EtOH	0,061	0,051	60
DY-845	847	160,000	876	EtOH	0,073	0,193	60
DY-865	863	190,000	896	EtOH	0,107	0,211	61

DY-350XL	347	17,000	608	PBS	0,477	0,401	63
DY-360XL	362	27,000	459	PBS	0,397	0,551	63
DY-370XL	368	17,000	475	PBS	0,460	0,733	64
DY-375XL	375	20,000	543	PBS	0,125	0,209	64
DY-376XL	387	25,400	549	PBS	0,507	0,575	65
DY-380XL	382	22,000	510	PBS	0,151	0,242	65
DY-395XL	397	22,000	572	PBS	0,210	0,541	66
DY-396XL	394	26,600	572	PBS	0,505	0,562	66
DY-475XL	494	90,000	512	EtOH	0,083	0,145	67
DY-480XL	504	50,000	631	EtOH	0,083	0,185	67
DY-481XL	521	50,000	649	EtOH	0,142	0,301	68
DY-485XL	488	50,000	559	EtOH	0,137	0,188	68
DY-494XL	497	50,000	639	EtOH	0,201	0,374	69
DY-510XL	512	50,000	590	EtOH	0,174	0,292	69
DY-511XL	510	50,000	595	PBS	0,230	0,309	70
DY-520XL	522	50,000	662	EtOH	0,217	0,425	70
DY-521XL	526	50,000	666	EtOH	0,184	0,379	71

Chemical structures in the supplement.

* The given values can change depending on the environment of the label (nature of conjugate and solvent).

Table with spectral characteristics of DY-Labels

label	absorption maximum* [nm]	molar absorbance* [l · mol ⁻¹ · cm ⁻¹]	emission maximum* [nm]	solvent	CF		page
					[280nm]	[260nm]	
DYQ-1	543	44,000	–	PBS	0,252	0,251	73
DYQ-2	641	105,000	–	PBS	0,126	0,183	73
DYQ-3	683	80,000	–	PBS	0,161	0,544	74
DYQ-4	766	180,000	–	PBS	n.d.	n.d.	74
DYQ-425	429	24,500	–	PBS	0,387	0,497	75
DYQ-505	506	44,000	–	PBS	0,277	0,281	75
DYQ-660	660	140,000	–	EtOH	n.d.	n.d.	76
DYQ-661	662	140,000	–	EtOH	0,080	0,052	76
DYQ-700	691	58,000	–	PBS	n.d.	n.d.	77
DY-485XLP4	478	50,000	557	EtOH	0,130	0,157	79
DY-520XLP4	514	50,000	665	EtOH	n.d.	n.d.	79
DY-547P4	557	150,000	571	PBS	0,083	0,046	80
DY-631P4	637	200,000	659	EtOH	0,075	0,065	80
DY-647P4	656	250,000	674	PBS	0,350	0,031	81
DY-679P4	684	200,000	703	PBS	0,108	0,097	81
DY-681P4	690	140,000	709	EtOH	0,093	0,129	82
DY-703P4	704	140,000	722	EtOH	0,153	0,225	82
DY-704P4	705	140,000	721	EtOH	0,129	0,163	83
DY-731P4	737	220,000	754	EtOH	0,047	0,051	83
DY-747P4	756	240,000	776	PBS	0,027	0,020	84
DY-781P4	784	170,000	795	EtOH	0,094	0,092	84
DY-800P4	783	250,000	795	PBS	0,041	0,040	85
FatRed	656	95,000	680	EtOH	n.d.	n.d.	87
Stain01	513	120,000	580	EtOH	n.d.	n.d.	87
DY-344IN	387	27,000	549	n.d.	n.d.	n.d.	88
V02-07027	477	70,000	522	EtOH	n.d.	n.d.	89
V02-08078	492	90,000	537	EtOH	n.d.	n.d.	89
V13-01184	481	80,000	526	EtOH	n.d.	n.d.	90
MitoDy-1	485	62,000	558	EtOH	n.d.	n.d.	90
DY-660-SL	665	20,000	–	EtOH	n.d.	n.d.	91
DY-840-S	840	65,000	–	EtOH	n.d.	n.d.	91
DY-660-X	662	80,000	693	Dichloromethane	n.d.	n.d.	92
DY-665-X	670	80,000	690	Dichloromethane	n.d.	n.d.	92


Chemical structures in the supplement.

* The given values can change depending on the environment of the label (nature of conjugate and solvent).

Aminoreactive Dyes

Primary amino groups are natural components of peptides and proteins which consist of amino acids. At least one amino group is positioned at the N-terminus of the peptide chain. Other biomolecules can be modified in such way that an amino group is incorporated without changing its activity or functionality. The most common examples are 5'-aminomodified DNA oligomers and cDNA bearing aminoallyl-dU-units.

There are several functional groups which react with primary amino groups forming a covalent bond between the dye and the amino containing (bio-)molecule under mild conditions not harmful to biologically active compounds. The most common reactive group being the NHS-ester (N-Hydroxysuccinimide-ester). At a pH between 8 and 9, the NHS-ester reacts specifically with primary amino groups. Depending on the purity (content of active NHS-ester), quantitative yields can be achieved.




Dyomics offers further aminoreactive groups too, e.g. TFP-ester (2,3,5,6-Tetrafluorophenyl-ester) or STP-ester (4-Sulfo-2,3,5,6-Tetrafluorophenyl-ester).

At Dyomics, special attention is paid to the activity of our labels. We carefully store and ship our DY-NHS-esters to supply our customers with high quality reagents. If the activity of the label suffered during shipment we replace the label for free.

Maleimides

As well as primary amino groups, thiol groups are preferred targets for fluorescent labels. There are various chemical moieties which form stable chemical bonds with thiol groups, including maleimides and iodoacetamides. The low influence of a neutral thiol-reactive label on the isoelectric point of a protein is an advantage over the labeling of amino groups in proteins and makes differential 2D-gel electrophoresis possible. Maleimides couple to thiol groups under very mild conditions at pH around 7.




Aminomodified labels

Labels with amino groups can be used for cross coupling to carboxy groups via a mild in-situ activation. Aminomodified labels can also be incorporated into a matrix during polymerization. Coupling reactions should be performed at pH 8.5 or higher since the dyes are delivered as ammonium salts.

Biotinylated labels and others



These modifications can be used to bind to streptavidin or avidin under physiological conditions, making it quite easy to detect the above proteins which are widely used in bioanalytics.

Dyomics offers a wide range of further reactive groups e.g. Hydrazide, Azide, Acetylene, Cadaverine. If you have need for a special modified label, please ask us.





Label | UV-excitation

DY-350

Absorption/emission max.:	553 nm / 432 nm (in PBS)		
Molar absorbance:	19,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	777.03	C ₁₉ H ₂₀ N ₂ O ₁₁ S ₂ * 2 C ₈ H ₂₀ N	350-00
NHS-ester	874.10	C ₂₃ H ₂₃ N ₃ O ₁₃ S ₂ * 2 C ₈ H ₂₀ N	350-01
Amino-derivative	689.86	C ₂₁ H ₂₇ N ₄ O ₁₀ S ₂ * C ₈ H ₂₀ N	350-02
Maleimide	899.15	C ₂₅ H ₂₆ N ₄ O ₁₂ S ₂ * 2 C ₈ H ₂₀ N	350-03
			



DY-351

Absorption/emission max.:	338 nm / 437 nm (in PBS)		
Molar absorbance:	33,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	454.43	C ₂₀ H ₁₉ N ₂ O ₇ S * Na	351-00
NHS-ester	551.50	C ₂₄ H ₂₂ N ₃ O ₉ S * Na	351-01
Amino-derivative	474.53	C ₂₂ H ₂₆ N ₄ O ₆ S	351-02
Maleimide	576.55	C ₂₆ H ₂₅ N ₄ O ₈ S * Na	351-03
Structure on request			





Label | Blue-green excitation



DY-405

Absorption/emission max.:	400 nm / 423 nm (in PBS)		
Molar absorbance:	32,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	695.59	C ₂₄ H ₂₀ NO ₁₃ S ₃ * 3 Na	405-00
NHS-ester	792.66	C ₂₈ H ₂₃ N ₂ O ₁₅ S ₃ * 3 Na	405-01
Amino-derivative	715.69	C ₂₆ H ₂₇ N ₃ O ₁₂ S ₃ * 2 Na	405-02
Maleimide	817.72	C ₃₀ H ₂₆ N ₃ O ₁₄ S ₃ * 3 Na	405-03
			

DY-410

Absorption/emission max.:	405 nm / 460 nm (in PBS)		
Molar absorbance:	34,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	455.78	C ₁₆ H ₁₅ ClNO ₉ S * Na	410-00
NHS-ester	552.87	C ₂₀ H ₁₈ ClN ₂ O ₁₁ S * Na	410-01
Amino-derivative	475.90	C ₁₈ H ₂₂ ClN ₃ O ₈ S	410-02
Maleimide	577.92	C ₂₂ H ₂₁ ClN ₃ O ₁₀ S * Na	410-03
			



DY-415

Absorption/emission max.:	418 nm / 463 nm (in Ethanol)		
Molar absorbance:	40,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	476.48	C ₂₀ H ₂₅ N ₂ O ₈ S * Na	415-00
NHS-ester	573.56	C ₂₄ H ₂₈ N ₃ O ₁₀ S * Na	415-01
Amino-derivative	496.59	C ₂₂ H ₃₂ N ₄ O ₇ S	415-02
Maleimide	598.61	C ₂₆ H ₃₁ N ₄ O ₉ S * Na	415-03
			



DY-430

Absorption/emission max.:	444 nm / 487 nm (PBS)		
Molar absorbance:	37,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	542.58	C ₂₅ H ₃₁ N ₂ O ₈ S * Na	430-00
NHS-ester	639.65	C ₂₉ H ₃₄ N ₃ O ₁₀ S * Na	430-01
Amino-derivative	562.68	C ₂₇ H ₃₈ N ₄ O ₇ S	430-02
Maleimide	664.70	C ₃₁ H ₃₇ N ₄ O ₉ S * Na	430-03



DY-431

Absorption/emission max.:	442 nm / 496 nm (in PBS)		
Molar absorbance:	35,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	644.63	C ₂₅ H ₃₀ N ₂ O ₁₁ S ₂ * 2 Na	431-00
NHS-ester	741.69	C ₂₉ H ₃₃ N ₃ O ₁₃ S ₂ * 2 Na	431-01
Amino-derivative	664.72	C ₂₇ H ₃₇ N ₄ O ₁₀ S ₂ * Na	431-02
Maleimide	766.75	C ₃₁ H ₃₆ N ₄ O ₁₂ S ₂ * 2 Na	431-03
			



DY-478

Absorption/emission max.:	478 nm / 518 nm (in Ethanol)		
Molar absorbance:	90,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	565.43	C ₂₇ H ₃₂ N ₃ O ₃ S * BF ₄	478-00
NHS-ester	662.50	C ₃₁ H ₃₅ N ₄ O ₅ S * BF ₄	478-01
Amino-derivative	592.62	C ₂₉ H ₃₉ N ₅ O ₂ S * 2 Cl	478-02
Maleimide	687.55	C ₃₃ H ₃₈ N ₅ O ₄ S * BF ₄	478-03
			



DY-490

Absorption/emission max.:	491 nm / 515 nm (in PBS)		
Molar absorbance:	73,000 M ⁻¹ cm ⁻¹		
Comments:			
– alternative to AlexaFluor™ 488			
– soluble in water, methanol			
– bright green, pH-stable emission			
– suitable for flow cytometry, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	677.61	C ₂₄ H ₂₀ N ₃ O ₁₃ S ₃ + Na	490-00
NHS-ester	1011.20	C ₂₈ H ₂₂ N ₄ O ₁₅ S ₃ * 2 C ₈ H ₂₀ N	490-01
Amino-derivative	925.78	C ₂₆ H ₂₉ N ₅ O ₁₂ S ₃ * 2 C ₂ F ₃ O ₂	490-02
Maleimide	799.74	C ₃₀ H ₂₆ N ₅ O ₁₄ S ₃ + Na	490-03
			


DY-495

Absorption/emission max.:	494 nm / 521 nm (in pH9-buffer)		
Molar absorbance:	78,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in methanol, ethanol, DMF, DMSO			
– bright solid state emission			
– suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	525.94	C ₂₇ H ₂₄ NO ₈ * Cl	495-00
NHS-ester	623.02	C ₃₁ H ₂₇ N ₂ O ₁₀ * Cl	495-01
Amino-derivative	604.49	C ₂₉ H ₃₁ N ₃ O ₇ * 2 Cl	495-02
Maleimide	725.62	C ₃₃ H ₃₀ N ₃ O ₉ * C ₂ F ₃ O ₂	495-03
			

DY-505



Absorption/emission max.:	507 nm / 528 nm (in Ethanol)		
Molar absorbance:	80,000 M ⁻¹ cm ⁻¹		
Comments:			
<ul style="list-style-type: none">– soluble in methanol, ethanol, DMF, DMSO– bright solid state emission– suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	523.97	C ₂₇ H ₂₆ N ₃ O ₆ · Cl	505-00
NHS-ester	621.05	C ₃₁ H ₂₉ N ₄ O ₈ · Cl	505-01
Amino-derivative	602.53	C ₂₉ H ₃₃ N ₅ O ₅ · 2 Cl	505-02
Maleimide	723.67	C ₃₃ H ₃₂ N ₅ O ₇ · C ₂ F ₃ O ₂	505-03
			


DYOMICS
Colours for Life





Label | Green excitation



DY-530

Absorption/emission max.:	533 nm / 554 nm (in PBS)		
Molar absorbance:	100,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	695.75	C ₃₁ H ₃₄ N ₃ O ₁₀ S ₂ * Na	530-00
NHS-ester	792.82	C ₃₅ H ₃₇ N ₄ O ₁₂ S ₂ * Na	530-01
Amino-derivative	715.85	C ₃₃ H ₄₁ N ₅ O ₉ S ₂	530-02
Maleimide	817.88	C ₃₇ H ₄₀ N ₅ O ₁₁ S ₂ * Na	530-03
			



DY-546

Absorption/emission max.:	549 nm / 569 nm (in PBS)		
Molar absorbance:	100,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	792.68	C ₃₃ H ₃₆ Cl ₂ N ₃ O ₁₀ S ₂ * Na	546-00
NHS-ester	889.75	C ₃₇ H ₃₉ Cl ₂ N ₄ O ₁₂ S ₂ * Na	546-01
Amino-derivative	812.78	C ₃₅ H ₄₃ Cl ₂ N ₅ O ₉ S ₂	546-02
Maleimide	914.80	C ₃₉ H ₄₂ Cl ₂ N ₅ O ₁₁ S ₂ * Na	546-03
			



DY-547P1

Absorption/emission max.:	551 nm / 565 nm (in PBS)		
Molar absorbance:	150,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	682.78	C ₃₂ H ₃₉ N ₂ O ₉ S ₂ * Na	547P1-00
NHS-ester	779.85	C ₃₆ H ₄₂ N ₃ O ₁₁ S ₂ * Na	547P1-01
Amino-derivative	702.88	C ₃₄ H ₄₆ N ₄ O ₈ S ₂	547P1-02
Maleimide	804.90	C ₃₈ H ₄₅ N ₄ O ₁₀ S ₂ * Na	547P1-03
			



DY-548P1

Absorption/emission max.:	554 nm / 568 nm (in PBS)		
Molar absorbance:	150,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	812.88	C ₃₄ H ₄₂ N ₂ O ₁₂ S ₃ * 2 Na	548P1-00
NHS-ester	909.95	C ₃₈ H ₄₅ N ₃ O ₁₄ S ₃ * 2 Na	548P1-01
Amino-derivative	832.98	C ₃₆ H ₄₉ N ₄ O ₁₁ S ₃ * Na	548P1-02
Maleimide	935.00	C ₄₀ H ₄₈ N ₄ O ₁₃ S ₃ * 2 Na	548P1-03
			



DY-549P1

Absorption/emission max.:	556 nm / 570 nm (in PBS)		
Molar absorbance:	150,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	942.99	C ₃₆ H ₄₅ N ₂ O ₁₅ S ₄ * 3 Na	549P1-00
NHS-ester	1040.06	C ₄₀ H ₄₈ N ₃ O ₁₇ S ₄ * 3 Na	549P1-01
Amino-derivative	963.09	C ₃₈ H ₅₂ N ₄ O ₁₄ S ₄ * 2 Na	549P1-02
Maleimide	1065.12	C ₄₂ H ₅₁ N ₄ O ₁₆ S ₄ * 3 Na	549P1-03
			



DY-550

Absorption/emission max.:	558 nm / 577 nm (in Ethanol)		
Molar absorbance:	120,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	667.76	C ₃₂ H ₄₂ N ₃ O ₉ S * Na	550-00
NHS-ester	764.83	C ₃₆ H ₄₅ N ₄ O ₁₁ S * Na	550-01
Amino-derivative	687.86	C ₃₄ H ₄₉ N ₅ O ₈ S	550-02
Maleimide	789.89	C ₃₈ H ₄₈ N ₅ O ₁₀ S * Na	550-03
			



DY-554

Absorption/emission max.:	544 nm / 570 nm (in Ethanol)		
Molar absorbance:	100,000 M ⁻¹ cm ⁻¹		
Comments:			
<ul style="list-style-type: none">– soluble in methanol, ethanol, DMF, DMSO– bright solid state emission– suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	580.08	C ₃₁ H ₃₄ N ₃ O ₆ * Cl	554-00
NHS-ester	677.16	C ₃₅ H ₃₇ N ₄ O ₈ * Cl	554-01
Amino-derivative	658.62	C ₃₃ H ₄₁ N ₅ O ₅ * 2 Cl	554-02
Maleimide	779.78	C ₃₇ H ₄₀ N ₅ O ₇ * C ₂ F ₃ O ₂	554-03
			



DY-555

Absorption/emission max.:	547 nm / 573 nm (in Ethanol)		
Molar absorbance:	100,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in methanol, ethanol, DMF, DMSO			
– bright solid state emission			
– suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	636.19	C ₃₅ H ₄₂ N ₃ O ₆ * Cl	555-00
NHS-ester	733.27	C ₃₉ H ₄₅ N ₄ O ₈ * Cl	555-01
Amino-derivative	714.72	C ₃₇ H ₄₉ N ₅ O ₅ * 2 Cl	555-02
Maleimide	721.86	C ₄₁ H ₄₇ N ₅ O ₇	555-03
			



DY-556

Absorption/emission max.:	548 nm / 574 nm (in Ethanol)		
Molar absorbance:	100,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in methanol, ethanol, DMF, DMSO			
– bright solid state emission			
– suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	715.80	C ₃₆ H ₄₂ N ₃ O ₉ S * Na	556-00
NHS-ester	812.88	C ₄₀ H ₄₅ N ₄ O ₁₁ S * Na	556-01
Amino-derivative	735.91	C ₃₈ H ₄₉ N ₅ O ₈ S	556-02
Maleimide	837.93	C ₄₂ H ₄₈ N ₅ O ₁₀ S * Na	556-03
			



DY-557

Absorption/emission max.:	556 nm / 576 nm (in Ethanol)		
Molar absorbance:	100,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	831.86	C ₃₇ H ₄₃ N ₃ O ₁₂ S ₂ * 2 Na	557-00
NHS-ester	928.93	C ₄₁ H ₄₆ N ₄ O ₁₄ S ₂ * 2 Na	557-01
Amino-derivative	851.96	C ₃₉ H ₅₀ N ₅ O ₁₁ S ₂ * Na	557-02
Maleimide	953.98	C ₄₃ H ₄₉ N ₅ O ₁₃ S ₂ * 2 Na	557-03
			



DY-560

Absorption/emission max.:	560 nm / 578 nm (in Ethanol)		
Molar absorbance:	120,000 M ⁻¹ cm ⁻¹		
Comments:			
<ul style="list-style-type: none">– soluble in methanol, ethanol, DMF, DMSO– bright solid state emission– suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	671.84	C ₃₃ H ₄₁ N ₃ O ₈ S ₂	560-00
NHS-ester	768.91	C ₃₇ H ₄₄ N ₄ O ₁₀ S ₂	560-01
Amino-derivative	750.38	C ₃₅ H ₄₈ N ₅ O ₇ S ₂ * Cl	560-02
Maleimide	793.97	C ₃₉ H ₄₇ N ₅ O ₉ S ₂	560-03
			


DY-580

Absorption/emission max.:	572 nm / 589 nm (in Ethanol)		
Molar absorbance:	100,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	571.06	C ₃₃ H ₃₁ N ₂ O ₅ * Cl	580-00
NHS-ester	745.70	C ₃₇ H ₃₄ N ₃ O ₇ * C ₂ F ₃ O ₂	580-01
Amino-derivative	649.60	C ₃₅ H ₃₈ N ₄ O ₄ * 2 Cl	580-02
Maleimide	770.75	C ₃₉ H ₃₇ N ₄ O ₆ * C ₂ F ₃ O ₂	580-03
			

DY-585



Absorption/emission max.:	584 nm / 604 nm (in Ethanol)		
Molar absorbance:	120,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	803.92	C ₃₉ H ₄₆ N ₃ O ₁₀ S ₂ * Na	585-00
NHS-ester	900.99	C ₄₃ H ₄₉ N ₄ O ₁₂ S ₂ * Na	585-01
Amino-derivative	824.02	C ₄₁ H ₅₃ N ₅ O ₉ S ₂	585-02
Maleimide	904.06	C ₄₅ H ₅₃ N ₅ O ₁₁ S ₂	585-03
			

DYOMICS
Colours for Life





Label | Orange excitation



DY-590

Absorption/emission max.:	581 nm / 600 nm (in Ethanol)		
Molar absorbance:	120,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, DMF, DMSO – bright solid state emission – suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	719.88	C ₃₇ H ₄₁ N ₃ O ₈ S ₂	590-00
NHS-ester	816.96	C ₄₁ H ₄₄ N ₄ O ₁₀ S ₂	590-01
Amino-derivative	798.42	C ₃₉ H ₄₈ N ₅ O ₇ S ₂ * Cl	590-02
Maleimide	842.01	C ₄₃ H ₄₇ N ₅ O ₉ S ₂	590-03
			



DY-591

Absorption/emission max.:	581 nm / 598 nm (in Ethanol)		
Molar absorbance:	120,000 M ⁻¹ cm ⁻¹		
Comments: – more hydrophilic than DY-590 – similar to Texas Red™			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	765.91	C ₃₈ H ₄₃ N ₃ O ₁₀ S ₂	591-00
NHS-ester	862.98	C ₄₂ H ₄₆ N ₄ O ₁₂ S ₂	591-01
Amino-derivative	844.45	C ₄₀ H ₅₀ N ₅ O ₉ S ₂ * Cl	591-02
Maleimide	888.04	C ₄₄ H ₄₉ N ₅ O ₁₁ S ₂	591-03
			



DY-594

Absorption/emission max.:	594 nm / 615 nm (in PBS)		
Molar absorbance:	92,000 M ⁻¹ cm ⁻¹		
Comments:	– similar to AlexaFluor™ 594		
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	937.06	C ₃₉ H ₄₄ N ₄ O ₁₅ S ₄	594-00
NHS-ester	1078.10	C ₄₃ H ₄₅ N ₅ O ₁₇ S ₄ * 2 Na	594-01
Amino-derivative	979.14	C ₄₁ H ₅₀ N ₆ O ₁₄ S ₄	594-02
Maleimide	1059.19	C ₄₅ H ₅₀ N ₆ O ₁₆ S ₄	594-03
			



DY-605


Absorption/emission max.:	600 nm / 624 nm (in PBS)		
Molar absorbance:	110,000 M ⁻¹ cm ⁻¹		
Comments:	– hydrophilic dye		
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	992.18	C ₄₃ H ₅₃ N ₅ O ₁₄ S ₄	605-00
NHS-ester	1111.22	C ₄₇ H ₅₅ N ₆ O ₁₆ S ₄ * Na	605-01
Amino-derivative	1034.25	C ₄₅ H ₅₉ N ₇ O ₁₃ S ₄	605-02
Maleimide	1136.27	C ₄₉ H ₅₈ N ₇ O ₁₅ S ₄ * Na	605-03
			

DY-610

Absorption/emission max.:	610 nm / 632 nm (in Ethanol)		
Molar absorbance:	80,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, DMF, DMSO – bright solid state emission – suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	667.76	C ₃₂ H ₄₂ N ₃ O ₉ S * Na	610-00
NHS-ester	764.82	C ₃₆ H ₄₅ N ₄ O ₁₁ S * Na	610-01
Amino-derivative	687.86	C ₃₄ H ₄₉ N ₅ O ₈ S	610-02
Maleimide	789.89	C ₃₈ H ₄₈ N ₅ O ₁₀ S * Na	610-03
			



DY-615

Absorption/emission max.:	623 nm / 643 nm (in Ethanol)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, ethanol, DMF, DMSO – bright solid state emission			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	578.73	C ₃₂ H ₃₈ N ₂ O ₆ S	615-00
NHS-ester	675.80	C ₃₆ H ₄₁ N ₃ O ₆ S	615-01
Amino-derivative	734.83	C ₃₄ H ₄₅ N ₄ O ₅ S * C ₂ F ₃ O ₂	615-02
Maleimide	700.86	C ₃₈ H ₄₄ N ₄ O ₇ S	615-03
			





Label | Red excitation



DY-630

Absorption/emission max.:	638 nm / 658 nm (in Ethanol)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, ethanol, DMF, DMSO – bright solid state emission – suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	634.84	C ₃₆ H ₄₆ N ₂ O ₆ S	630-00
NHS-ester	731.92	C ₄₀ H ₄₉ N ₃ O ₈ S	630-01
Amino-derivative	790.93	C ₃₈ H ₅₃ N ₄ O ₅ S * C ₂ F ₃ O ₂	630-02
Maleimide	756.97	C ₄₂ H ₅₂ N ₄ O ₇ S	630-03
			



DY-631

Absorption/emission max.:	637 nm / 657 nm (in Ethanol)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO – bright solid state emission – suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	736.88	C ₃₆ H ₄₅ N ₂ O ₉ S ₂ * Na	631-00
NHS-ester	833.95	C ₄₀ H ₄₈ N ₃ O ₁₁ S ₂ * Na	631-01
Amino-derivative	756.99	C ₃₈ H ₅₂ N ₄ O ₈ S ₂	631-02
Maleimide	859.01	C ₄₂ H ₅₁ N ₄ O ₁₀ S ₂ * Na	631-03
			



DY-632

Absorption/emission max.:	636 nm / 658 nm (in Ethanol)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments: <ul style="list-style-type: none">– soluble in methanol, ethanol, DMF, DMSO– bright solid state emission– suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	852.96	C ₃₇ H ₄₆ N ₂ O ₁₂ S ₃ * 2 Na	632-00
NHS-ester	950.03	C ₄₁ H ₄₉ N ₃ O ₁₄ S ₃ * 2 Na	632-01
Amino-derivative	873.05	C ₃₉ H ₅₃ N ₄ O ₁₁ S ₃ * Na	632-02
Maleimide	975.08	C ₄₃ H ₅₂ N ₄ O ₁₃ S ₃ * 2 Na	632-03
			



DY-633

Absorption/emission max.:	638 nm / 658 nm (in Ethanol)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in water, methanol, DMF, DMSO			
– bright solid state emission			
– suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	750.91	C ₃₇ H ₄₇ N ₂ O ₉ S ₂ * Na	633-00
NHS-ester	847.98	C ₄₁ H ₅₀ N ₃ O ₁₁ S ₂ * Na	633-01
Amino-derivative	771.02	C ₃₉ H ₅₄ N ₄ O ₈ S ₂	633-02
Maleimide	873.04	C ₄₃ H ₅₃ N ₄ O ₁₀ S ₂ * Na	633-03
			



DY-634

Absorption/emission max.:	636 nm / 657 nm (in Ethanol)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in methanol, ethanol, DMF, DMSO			
– bright solid state emission			
– suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	969.03	C ₃₈ H ₄₇ N ₂ O ₁₅ S ₄ * 3 Na	634-00
NHS-ester	1066.10	C ₄₂ H ₅₀ N ₃ O ₁₇ S ₄ * 3 Na	634-01
Amino-derivative	989.12	C ₄₀ H ₅₄ N ₄ O ₁₄ S ₄ * 2 Na	634-02
Maleimide	1091.16	C ₄₄ H ₅₃ N ₄ O ₁₆ S ₄ * 3 Na	634-03
			



DY-635

Absorption/emission max.:	648 nm / 670 nm (in Ethanol)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments:			
<ul style="list-style-type: none">– soluble in methanol, ethanol, DMF, DMSO– bright solid state emission– suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	658.86	C ₃₈ H ₄₆ N ₂ O ₆ S	635-00
NHS-ester	755.93	C ₄₂ H ₄₉ N ₃ O ₆ S	635-01
Amino-derivative	737.41	C ₄₀ H ₅₃ N ₄ O ₅ S * Cl	635-02
Maleimide	780.99	C ₄₄ H ₅₂ N ₄ O ₇ S	635-03
			



DY-636

Absorption/emission max.:	647 nm / 670 nm (in Ethanol)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO – bright solid state emission – suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	760.91	C ₃₈ H ₄₅ N ₂ O ₉ S ₂ * Na	636-00
NHS-ester	857.98	C ₄₂ H ₄₈ N ₃ O ₁₁ S ₂ * Na	636-01
Amino-derivative	781.01	C ₄₀ H ₅₂ N ₄ O ₈ S ₂	636-02
Maleimide	883.04	C ₄₄ H ₅₁ N ₄ O ₁₀ S ₂ * Na	636-03
			



DY-641

Absorption/emission max.:	624 nm / 654 nm (in Ethanol)		
Molar absorbance:	80,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	803.81	C ₃₅ H ₃₉ N ₃ O ₁₂ S ₂ * 2 Na	641-00
NHS-ester	900.88	C ₃₉ H ₄₂ N ₄ O ₁₄ S ₂ * 2 Na	641-01
Amino-derivative	823.91	C ₃₇ H ₄₆ N ₅ O ₁₁ S ₂ * Na	641-02
Maleimide	925.93	C ₄₁ H ₄₅ N ₅ O ₁₃ S ₂ * 2 Na	641-03
			



DY-643

Absorption/emission max.:	626 nm / 646 nm (in Ethanol)		
Molar absorbance:	80,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	833.83	C ₃₆ H ₄₁ N ₃ O ₁₃ S ₂ * 2 Na	643-00
NHS-ester	930.90	C ₄₀ H ₄₄ N ₄ O ₁₅ S ₂ * 2 Na	643-01
Amino-derivative	853.93	C ₃₈ H ₄₈ N ₅ O ₁₂ S ₂ * Na	643-02
Maleimide	955.96	C ₄₂ H ₄₇ N ₅ O ₁₄ S ₂ * 2 Na	643-03
			



DY-647P1

Absorption/emission max.:	648 nm / 666 nm (in PBS)		
Molar absorbance:	250,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	708.82	C ₃₄ H ₄₁ N ₂ O ₉ S ₂ * Na	647P1-00
NHS-ester	805.89	C ₃₈ H ₄₄ N ₃ O ₁₁ S ₂ * Na	647P1-01
Amino-derivative	728.92	C ₃₆ H ₄₈ N ₄ O ₈ S ₂	647P1-02
Maleimide	830.94	C ₄₀ H ₄₇ N ₄ O ₁₀ S ₂ * Na	647P1-03
			



DY-648P1

Absorption/emission max.:	651 nm / 669 nm (in PBS)		
Molar absorbance:	250,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	838.92	C ₃₆ H ₄₄ N ₂ O ₁₂ S ₃ * 2 Na	648P1-00
NHS-ester	935.99	C ₄₀ H ₄₇ N ₃ O ₁₄ S ₃ * 2 Na	648P1-01
Amino-derivative	859.02	C ₃₈ H ₅₁ N ₄ O ₁₁ S ₃ * Na	648P1-02
Maleimide	961.04	C ₄₂ H ₅₀ N ₄ O ₁₃ S ₃ * 2 Na	648P1-03
			



DY-649P1

Absorption/emission max.:	654 nm / 672 nm (in PBS)		
Molar absorbance:	250,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol – very hydrophilic – three negative charges – usually combined with DY-549P1			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	969.03	C ₃₈ H ₄₇ N ₂ O ₁₅ S ₄ * 3 Na	649P1-00
NHS-ester	1066.10	C ₄₂ H ₅₀ N ₃ O ₁₇ S ₄ * 3 Na	649P1-01
Amino-derivative	989.13	C ₄₀ H ₅₄ N ₄ O ₁₄ S ₄ * 2 Na	649P1-02
Maleimide	1091.16	C ₄₄ H ₅₃ N ₄ O ₁₆ S ₄ * 3 Na	649P1-03
			



DY-650

Absorption/emission max.:	656 nm / 676 nm (in Ethanol)		
Molar absorbance:	220,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, ethanol, DMF, DMSO – bright solid state emission – suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	686.92	C ₄₀ H ₅₀ N ₂ O ₆ S	650-00
NHS-ester	783.99	C ₄₄ H ₅₃ N ₃ O ₈ S	650-01
Amino-derivative	765.46	C ₄₂ H ₅₇ N ₄ O ₅ S * Cl	650-02
Maleimide	809.05	C ₄₆ H ₅₆ N ₄ O ₇ S	650-03
			



DY-651

Absorption/emission max.:	655 nm / 677 nm (in Ethanol)		
Molar absorbance:	220,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO – bright solid state emission – suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	788.96	C ₄₀ H ₄₉ N ₂ O ₉ S ₂ * Na	651-00
NHS-ester	886.04	C ₄₄ H ₅₂ N ₃ O ₁₁ S ₂ * Na	651-01
Amino-derivative	809.06	C ₄₂ H ₅₆ N ₄ O ₈ S ₂	651-02
Maleimide	911.09	C ₄₆ H ₅₅ N ₄ O ₁₀ S ₂ * Na	651-03
			



DY-652

Absorption/emission max.:	653 nm / 676 nm (in Ethanol)		
Molar absorbance:	220,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	905.03	C ₄₁ H ₅₀ N ₂ O ₁₂ S ₃ * 2 Na	652-00
NHS-ester	1002.11	C ₄₅ H ₅₃ N ₃ O ₁₄ S ₃ * 2 Na	652-01
Amino-derivative	925.13	C ₄₃ H ₅₇ N ₄ O ₁₁ S ₃ * Na	652-02
Maleimide	1027.16	C ₄₇ H ₅₆ N ₄ O ₁₃ S ₃ * 2 Na	652-03
			

DY-654

Absorption/emission max.:	653 nm / 677 nm (in Ethanol)		
Molar absorbance:	220,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1007.08	C ₄₁ H ₄₉ N ₂ O ₁₅ S ₄ * 3 Na	654-00
NHS-ester	1104.13	C ₄₅ H ₅₂ N ₃ O ₁₇ S ₄ * 3 Na	654-01
Amino-derivative	1027.18	C ₄₃ H ₅₆ N ₄ O ₁₄ S ₄ * 2 Na	654-02
Maleimide	1129.21	C ₄₇ H ₅₅ N ₄ O ₁₆ S ₄ * 3 Na	654-03
			

DY-660P1



Absorption/emission max.:	663 nm / 682 nm (in PBS)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	991.02	C ₄₀ H ₄₅ N ₂ O ₁₅ S ₄ * 3 Na	660P1-00
NHS-ester	1088.09	C ₄₄ H ₄₈ N ₃ O ₁₇ S ₄ * 3 Na	660P1-01
Amino-derivative	1011.12	C ₄₂ H ₅₂ N ₄ O ₁₄ S ₄ * 2 Na	660P1-02
Maleimide	1113.14	C ₄₆ H ₅₁ N ₄ O ₁₆ S ₄ * 3 Na	660P1-03
			

DYOMICS
Colours for Life





Label | Far red excitation



DY-675

Absorption/emission max.:	675 nm / 699 nm (in Ethanol)		
Molar absorbance:	180,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, ethanol, DMF, DMSO – bright solid state emission – suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	706.91	C ₄₂ H ₄₆ N ₂ O ₆ S	675-00
NHS-ester	803.98	C ₄₆ H ₄₉ N ₃ O ₆ S	675-01
Amino-derivative	785.45	C ₄₄ H ₅₃ N ₄ O ₅ S * Cl	675-02
Maleimide	829.04	C ₄₈ H ₅₂ N ₄ O ₇ S	675-03
			



DY-676

Absorption/emission max.:	675 nm / 699 nm (in Ethanol)		
Molar absorbance:	180,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO – bright solid state emission – suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	808.95	C ₄₂ H ₄₅ N ₂ O ₉ S ₂ * Na	676-00
NHS-ester	906.02	C ₄₆ H ₄₈ N ₃ O ₁₁ S ₂ * Na	676-01
Amino-derivative	829.06	C ₄₄ H ₅₂ N ₄ O ₈ S ₂	676-02
Maleimide	931.08	C ₄₈ H ₅₁ N ₄ O ₁₀ S ₂ * Na	676-03
			

DY-677

Absorption/emission max.:	674 nm / 698 nm (in Ethanol)		
Molar absorbance:	180,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	925.02	C ₄₃ H ₄₆ N ₂ O ₁₂ S ₃ * 2 Na	677-00
NHS-ester	1022.10	C ₄₇ H ₄₉ N ₃ O ₁₄ S ₃ * 2 Na	677-01
Amino-derivative	945.12	C ₄₅ H ₅₃ N ₄ O ₁₁ S ₃ * Na	677-02
Maleimide	1047.15	C ₄₉ H ₅₂ N ₄ O ₁₃ S ₃ * 2 Na	677-03
			



DY-678

Absorption/emission max.:	674 nm / 694 nm (in Ethanol)		
Molar absorbance:	180,000 M ⁻¹ cm ⁻¹		
Comments: – very hydrophilic			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1027.07	C ₄₃ H ₄₅ N ₂ O ₁₅ S ₄ * 3 Na	678-00
NHS-ester	1124.14	C ₄₇ H ₄₈ N ₃ O ₁₇ S ₄ * 3 Na	678-01
Amino-derivative	1047.17	C ₄₅ H ₅₂ N ₄ O ₁₄ S ₄ * 2 Na	678-02
Maleimide	1149.18	C ₄₉ H ₅₁ N ₄ O ₁₆ S ₄ * 3 Na	678-03
			



DY-679P1

Absorption/emission max.:	679 nm / 697 nm (in PBS)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF-water-mixtures			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1143.14	C ₄₄ H ₄₆ N ₂ O ₁₈ S ₅ * 4 Na	679P1-00
NHS-ester	1240.21	C ₄₈ H ₄₉ N ₃ O ₂₀ S ₅ * 4 Na	679P1-01
Amino-derivative	1163.24	C ₄₆ H ₅₃ N ₄ O ₁₇ S ₅ * 3 Na	679P1-02
Maleimide	1265.27	C ₅₀ H ₅₂ N ₄ O ₁₉ S ₅ * 4 Na	679P1-03



DY-680

Absorption/emission max.:	691 nm / 709 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, ethanol, DMF, DMSO – bright solid state emission – suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	634.84	C ₃₆ H ₄₆ N ₂ O ₆ S	680-00
NHS-ester	731.92	C ₄₀ H ₄₉ N ₃ O ₆ S	680-01
Amino-derivative	713.38	C ₃₈ H ₅₃ N ₄ O ₅ S * Cl	680-02
Maleimide	756.97	C ₄₂ H ₅₂ N ₄ O ₇ S	680-03
			



DY-681

Absorption/emission max.:	692 nm / 709 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in water, methanol, DMF, DMSO			
– bright solid state emission			
– suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	736.88	C ₃₆ H ₄₅ N ₂ O ₉ S ₂ * Na	681-00
NHS-ester	833.95	C ₄₀ H ₄₈ N ₃ O ₁₁ S ₂ * Na	681-01
Amino-derivative	756.99	C ₃₈ H ₅₂ N ₄ O ₈ S ₂	681-02
Maleimide	859.01	C ₄₂ H ₅₁ N ₄ O ₁₀ S ₂ * Na	681-03
			


DY-682

Absorption/emission max.:	692 nm / 709 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: - soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	852.96	C ₃₇ H ₄₆ N ₂ O ₁₂ S ₃ * 2 Na	682-00
NHS-ester	950.03	C ₄₁ H ₄₉ N ₃ O ₁₄ S ₃ * 2 Na	682-01
Amino-derivative	873.06	C ₃₉ H ₅₃ N ₄ O ₁₁ S ₃ * Na	682-02
Maleimide	975.08	C ₄₃ H ₅₂ N ₄ O ₁₃ S ₃ * 2 Na	682-03
			

DY-684



Absorption/emission max.:	685 nm / 699 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1065.10	C ₄₆ H ₄₇ N ₂ O ₁₅ S ₄ * 3 Na	684-00
NHS-ester	1162.17	C ₅₀ H ₅₀ N ₃ O ₁₇ S ₄ * 3 Na	684-01
Amino-derivative	1085.20	C ₄₈ H ₅₃ N ₄ O ₁₄ S ₄ * 2 Na	684-02
Maleimide	1187.22	C ₅₂ H ₅₂ N ₄ O ₁₆ S ₄ * 3 Na	684-03
			

DYOMICS
Colours for Life





 **Label** | NIR-excitation



DY-700

Absorption/emission max.:	707 nm / 728 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, ethanol, DMF, DMSO – bright solid state emission – suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	668.86	C ₃₉ H ₄₄ N ₂ O ₆ S	700-00
NHS-ester	765.93	C ₄₃ H ₄₇ N ₃ O ₈ S	700-01
Amino-derivative	747.40	C ₄₁ H ₅₁ N ₄ O ₅ S * Cl	700-02
Maleimide	790.99	C ₄₅ H ₅₀ N ₄ O ₇ S	700-03
			



DY-701

Absorption/emission max.:	709 nm / 730 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in water, methanol, DMF, DMSO			
– bright solid state emission			
– suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	770.90	C ₃₉ H ₄₃ N ₂ O ₉ S ₂ * Na	701-00
NHS-ester	867.97	C ₄₃ H ₄₆ N ₃ O ₁₁ S ₂ * Na	701-01
Amino-derivative	791.01	C ₄₁ H ₅₀ N ₄ O ₈ S ₂	701-02
Maleimide	893.03	C ₄₅ H ₄₉ N ₄ O ₁₀ S ₂ * Na	701-03
			



DY-702

Absorption/emission max.:	709 nm / 728 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	886.96	C ₄₀ H ₄₄ N ₂ O ₁₂ S ₃ * 2 Na	702-00
NHS-ester	984.03	C ₄₄ H ₄₇ N ₃ O ₁₄ S ₃ * 2 Na	702-01
Amino-derivative	907.06	C ₄₂ H ₅₁ N ₄ O ₁₁ S ₃ * Na	702-02
Maleimide	1009.08	C ₄₆ H ₅₀ N ₄ O ₁₃ S ₃ * 2 Na	702-03
			



DY-703

Absorption/emission max.:	705 nm / 721 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	800.93	C ₄₀ H ₄₅ N ₂ O ₁₀ S ₂ * Na	703-00
NHS-ester	898.00	C ₄₄ H ₄₈ N ₃ O ₁₂ S ₂ * Na	703-01
Amino-derivative	821.03	C ₄₂ H ₅₂ N ₄ O ₉ S ₂	703-02
Maleimide	923.06	C ₄₆ H ₅₁ N ₄ O ₁₁ S ₂ * Na	703-03
			



DY-704

Absorption/emission max.:	706 nm / 721 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – very hydrophilic			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	902.97	C ₄₀ H ₄₄ N ₂ O ₁₃ S ₃ * 2 Na	704-00
NHS-ester	1000.05	C ₄₄ H ₄₇ N ₃ O ₁₅ S ₃ * 2 Na	704-01
Amino-derivative	923.07	C ₄₂ H ₅₁ N ₄ O ₁₂ S ₃ * Na	704-02
Maleimide	1025.10	C ₄₆ H ₅₀ N ₄ O ₁₄ S ₃ * 2 Na	704-03
			



DY-705

Absorption/emission max.:	710 nm / 732 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	770.89	C ₃₉ H ₄₃ N ₂ O ₉ S ₂ * Na	705-00
NHS-ester	867.96	C ₄₃ H ₄₆ N ₃ O ₁₁ S ₂ * Na	705-01
Amino-derivative	790.99	C ₄₁ H ₅₀ N ₄ O ₈ S ₂	705-02
Maleimide	893.01	C ₄₅ H ₄₉ N ₄ O ₁₀ S ₂ * Na	705-03
			



DY-706

Absorption/emission max.:	710 nm / 733 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	872.93	C ₃₉ H ₄₂ N ₂ O ₁₂ S ₃ * 2 Na	706-00
NHS-ester	970.00	C ₄₃ H ₄₅ N ₃ O ₁₄ S ₃ * 2 Na	706-01
Amino-derivative	893.03	C ₄₁ H ₄₉ N ₄ O ₁₁ S ₃ * Na	706-02
Maleimide	995.06	C ₄₅ H ₄₈ N ₄ O ₁₃ S ₃ * 2 Na	706-03
			



DY-720

Absorption/emission max.:	726 nm / 744 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	882.03	C ₄₅ H ₅₂ N ₃ O ₁₀ S ₂ * Na	720-00
NHS-ester	979.10	C ₄₉ H ₅₅ N ₄ O ₁₂ S ₂ * Na	720-01
Amino-derivative	902.13	C ₄₇ H ₅₉ N ₅ O ₉ S ₂	720-02
Maleimide	1004.15	C ₅₁ H ₅₈ N ₅ O ₁₁ S ₂ * Na	720-03
			



DY-730

Absorption/emission max.:	734 nm / 755 nm (in Ethanol)		
Molar absorbance:	240,000 M ⁻¹ cm ⁻¹		
Comments:			
<div><div></div><div>– soluble in methanol, ethanol, DMF, DMSO</div><div>– bright solid state emission</div><div>– suitable for microarray experiments, FisH microscopy, gel electrophoresis</div></div>			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	660.88	C ₃₈ H ₄₈ N ₂ O ₆ S	730-00
NHS-ester	757.96	C ₄₂ H ₅₁ N ₃ O ₈ S	730-01
Amino-derivative	739.42	C ₄₀ H ₅₅ N ₄ O ₅ S * Cl	730-02
Maleimide	783.01	C ₄₄ H ₅₄ N ₄ O ₇ S	730-03
<div></div>		<div></div>	



DY-731

Absorption/emission max.:	736 nm / 755 nm (in Ethanol)		
Molar absorbance:	240,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in water, methanol, DMF, DMSO			
– bright solid state emission			
– suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	762.92	C ₃₈ H ₄₇ N ₂ O ₉ S ₂ * Na	731-00
NHS-ester	859.99	C ₄₂ H ₅₀ N ₃ O ₁₁ S ₂ * Na	731-01
Amino-derivative	783.03	C ₄₀ H ₅₄ N ₄ O ₈ S ₂	731-02
Maleimide	885.05	C ₄₄ H ₅₃ N ₄ O ₁₀ S ₂ * Na	731-03
			



DY-732

Absorption/emission max.:	735 nm / 756 nm (in Ethanol)		
Molar absorbance:	240,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO – bright solid state emission – suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	879.00	C ₃₉ H ₄₈ N ₂ O ₁₂ S ₃ * 2 Na	732-00
NHS-ester	976.07	C ₄₃ H ₅₁ N ₃ O ₁₄ S ₃ * 2 Na	732-01
Amino-derivative	899.10	C ₄₁ H ₅₅ N ₄ O ₁₁ S ₃ * Na	732-02
Maleimide	1001.12	C ₄₅ H ₅₄ N ₄ O ₁₃ S ₃ * 2 Na	732-03
			



DY-734

Absorption/emission max.:	733 nm / 755 nm (in Ethanol)		
Molar absorbance:	240,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	995.07	C ₄₀ H ₄₉ N ₂ O ₁₅ S ₄ * 3 Na	734-00
NHS-ester	1092.14	C ₄₄ H ₅₂ N ₃ O ₁₇ S ₄ * 3 Na	734-01
Amino-derivative	1015.19	C ₄₂ H ₅₆ N ₄ O ₁₄ S ₄ * 2 Na	734-02
Maleimide	1117.19	C ₄₆ H ₅₅ N ₄ O ₁₆ S ₄ * 3 Na	734-03
			



DY-736

Absorption/emission max.:	743 nm / 764 nm (in Ethanol)		
Molar absorbance:	240,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	786.93	C ₄₀ H ₄₇ N ₂ O ₉ S ₂ * Na	736-00
NHS-ester	884.00	C ₄₄ H ₅₀ N ₃ O ₁₁ S ₂ * Na	736-01
Amino-derivative	807.03	C ₄₂ H ₅₄ N ₄ O ₈ S ₂	736-02
Maleimide	909.05	C ₄₆ H ₅₃ N ₄ O ₁₀ S ₂ * Na	736-03
			



DY-747P1

Absorption/emission max.:	747 nm / 769 nm (in PBS)		
Molar absorbance:	240,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	734.85	C ₃₆ H ₄₃ N ₂ O ₉ S ₂ * Na	747P1-00
NHS-ester	831.93	C ₄₀ H ₄₆ N ₃ O ₁₁ S ₂ * Na	747P1-01
Amino-derivative	754.96	C ₃₈ H ₅₀ N ₄ O ₈ S ₂	747P1-02
Maleimide	856.98	C ₄₂ H ₄₉ N ₄ O ₁₀ S ₂ * Na	747P1-03
			



DY-749P1

Absorption/emission max.:	747 nm / 769 nm (in PBS)		
Molar absorbance:	240,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	995.05	C ₄₀ H ₄₉ N ₂ O ₁₅ S ₄ * 3 Na	749P1-00
NHS-ester	1092.12	C ₄₄ H ₅₂ N ₃ O ₁₇ S ₄ * 3 Na	749P1-01
Amino-derivative	1015.15	C ₄₂ H ₅₆ N ₄ O ₁₄ S ₄ * 2 Na	749P1-02
Maleimide	1117.18	C ₄₆ H ₅₅ N ₄ O ₁₆ S ₄ * 3 Na	749P1-03
			



DY-750

Absorption/emission max.:	751 nm / 774 nm (in Ethanol)		
Molar absorbance:	270,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in methanol, ethanol, DMF, DMSO			
– bright solid state emission			
– suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	712.96	C ₄₂ H ₅₂ N ₂ O ₆ S	750-00
NHS-ester	810.03	C ₄₆ H ₅₅ N ₃ O ₈ S	750-01
Amino-derivative	791.49	C ₄₄ H ₅₉ N ₄ O ₅ S * Cl	750-02
Maleimide	835.08	C ₄₈ H ₅₈ N ₄ O ₇ S	750-03
			



DY-751

Absorption/emission max.:	752 nm / 772 nm (in Ethanol)		
Molar absorbance:	270,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO – bright solid state emission – suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	814.99	C ₄₂ H ₅₁ N ₂ O ₉ S ₂ * Na	751-00
NHS-ester	912.07	C ₄₆ H ₅₄ N ₃ O ₁₁ S ₂ * Na	751-01
Amino-derivative	835.10	C ₄₄ H ₅₈ N ₄ O ₈ S ₂	751-02
Maleimide	937.13	C ₄₈ H ₅₇ N ₄ O ₁₀ S ₂ * Na	751-03
			



DY-752

Absorption/emission max.:	750 nm / 771 nm (in Ethanol)		
Molar absorbance:	270,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	931.07	C ₄₃ H ₅₂ N ₂ O ₁₂ S ₃ * 2 Na	752-00
NHS-ester	1028.15	C ₄₇ H ₅₅ N ₃ O ₁₄ S ₃ * 2 Na	752-01
Amino-derivative	951.17	C ₄₅ H ₅₉ N ₄ O ₁₁ S ₃ * Na	752-02
Maleimide	1053.20	C ₄₉ H ₅₈ N ₄ O ₁₃ S ₃ * 2 Na	752-03
			



DY-754

Absorption/emission max.:	748 nm / 771 nm (in Ethanol)		
Molar absorbance:	270,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1033.12	C ₄₃ H ₅₁ N ₂ O ₁₅ S ₄ * 3 Na	754-00
NHS-ester	1130.19	C ₄₇ H ₅₄ N ₃ O ₁₇ S ₄ * 3 Na	754-01
Amino-derivative	1053.22	C ₄₅ H ₅₈ N ₄ O ₁₄ S ₄ * 2 Na	754-02
Maleimide	1155.25	C ₄₉ H ₅₇ N ₄ O ₁₆ S ₄ * 3 Na	754-03
			



DY-765

Absorption/emission max.:	765 nm / 785 nm (in Ethanol)		
Molar absorbance:	130,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1058.15	C ₅₁ H ₅₃ N ₃ O ₁₃ S ₃ * 2 Na	765-00
NHS-ester	1155.23	C ₅₅ H ₅₆ N ₄ O ₁₅ S ₃ * 2 Na	765-01
Amino-derivative	1078.25	C ₅₃ H ₆₀ N ₅ O ₁₂ S ₃ * Na	765-02
Maleimide	1180.28	C ₅₇ H ₅₉ N ₅ O ₁₄ S ₃ * 2 Na	765-03
			



DY-776

Absorption/emission max.:	772 nm / 787 nm (in Ethanol)		
Molar absorbance:	240,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO – bright solid state emission – suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	834.98	C ₄₄ H ₄₇ N ₂ O ₉ S ₂ * Na	776-00
NHS-ester	932.07	C ₄₈ H ₅₀ N ₃ O ₁₁ S ₂ * Na	776-01
Amino-derivative	855.09	C ₄₆ H ₅₄ N ₄ O ₈ S ₂	776-02
Maleimide	957.12	C ₅₀ H ₅₃ N ₄ O ₁₀ S ₂ * Na	776-03
			



DY-777

Absorption/emission max.:	770 nm / 788 nm (in Ethanol)		
Molar absorbance:	240,000 M ⁻¹ cm ⁻¹		
Comments: – hydrophilic			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	951.06	C ₄₅ H ₄₈ N ₂ O ₁₂ S ₃ * 2 Na	777-00
NHS-ester	1048.14	C ₄₉ H ₅₁ N ₃ O ₁₄ S ₃ * 2 Na	777-01
Amino-derivative	971.16	C ₄₇ H ₅₅ N ₄ O ₁₁ S ₃ * Na	777-02
Maleimide	1073.19	C ₅₁ H ₅₄ N ₄ O ₁₃ S ₃ * 2 Na	777-03
			



DY-778

Absorption/emission max.:	767 nm / 787 nm (in Ethanol)		
Molar absorbance:	240,000 M ⁻¹ cm ⁻¹		
Comments: – very hydrophilic			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1053.11	C ₄₅ H ₄₇ N ₂ O ₁₅ S ₄ * 3 Na	778-00
NHS-ester	1150.18	C ₄₉ H ₅₀ N ₃ O ₁₇ S ₄ * 3 Na	778-01
Amino-derivative	1073.21	C ₄₇ H ₅₄ N ₄ O ₁₄ S ₄ * 2 Na	778-02
Maleimide	1175.24	C ₅₁ H ₅₃ N ₄ O ₁₆ S ₄ * 3 Na	778-03
			



DY-780

Absorption/emission max.:	783 nm / 799 nm (in Ethanol)		
Molar absorbance:	170,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in methanol, ethanol, DMF, DMSO			
– bright solid state emission			
– suitable for microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	660.88	C ₃₈ H ₄₈ N ₂ O ₆ S	780-00
NHS-ester	757.95	C ₄₂ H ₅₁ N ₃ O ₆ S	780-01
Amino-derivative	816.97	C ₃₁ H ₅₅ N ₄ O ₅ S * C ₂ F ₃ O ₂	780-02
Maleimide	783.01	C ₄₄ H ₅₄ N ₄ O ₇ S	780-03
			



DY-781

Absorption/emission max.:	784 nm / 796 nm (in Ethanol)		
Molar absorbance:	170,000 M ⁻¹ cm ⁻¹		
Comments:			
<ul style="list-style-type: none">– soluble in water, methanol, DMF, DMSO– bright solid state emission– suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	762.92	C ₃₈ H ₄₇ N ₂ O ₉ S ₂ * Na	781-00
NHS-ester	860.00	C ₄₂ H ₅₀ N ₃ O ₁₁ S ₂ * Na	781-01
Amino-derivative	783.03	C ₄₀ H ₅₄ N ₄ O ₈ S ₂	781-02
Maleimide	885.05	C ₄₄ H ₅₃ N ₄ O ₁₀ S ₂ * Na	781-03
			



DY-782

Absorption/emission max.:	785 nm / 794 nm (in Ethanol)		
Molar absorbance:	170,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in water, methanol			
– two negative charges			
– enhanced water solubility and polarity			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	879.00	C ₃₉ H ₄₈ N ₂ O ₁₂ S ₃ * 2 Na	782-00
NHS-ester	976.07	C ₄₃ H ₅₁ N ₃ O ₁₄ S ₃ * 2 Na	782-01
Amino-derivative	899.10	C ₄₁ H ₅₅ N ₄ O ₁₁ S ₃ * Na	782-02
Maleimide	1001.12	C ₄₅ H ₅₄ N ₄ O ₁₃ S ₃ * 2 Na	782-03
			



DY-784

Absorption/emission max.:	779 nm / 780 nm (in Ethanol)		
Molar absorbance:	170,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1091.14	C ₄₈ H ₄₉ N ₂ O ₁₅ S ₄ * 3 Na	784-00
NHS-ester	1188.21	C ₅₂ H ₅₂ N ₃ O ₁₇ S ₄ * 3 Na	784-01
Amino-derivative	1111.24	C ₅₀ H ₅₆ N ₄ O ₁₄ S ₄ * 2 Na	784-02
Maleimide	1213.26	C ₅₄ H ₅₅ N ₄ O ₁₆ S ₄ * 3 Na	784-03
			



DY-800

Absorption/emission max.:	777 nm / 791 nm (in Ethanol)		
Molar absorbance:	280,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF – hydrophilic			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	953.08	C ₄₅ H ₅₀ N ₂ O ₁₂ S ₃ * 2 Na	800-00
NHS-ester	1050.15	C ₄₉ H ₅₃ N ₃ O ₁₄ S ₃ * 2 Na	800-01
Amino-derivative	973.18	C ₄₇ H ₅₇ N ₄ O ₁₁ S ₃ * Na	800-02
Maleimide	1075.21	C ₅₁ H ₅₆ N ₄ O ₁₃ S ₃ * 2 Na	800-03
			



DY-805

Absorption/emission max.:	800 nm / 824 nm (in Ethanol)		
Molar absorbance:	150,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	866.97	C ₄₄ H ₄₇ N ₂ O ₁₁ S ₂ * Na	805-00
NHS-ester	964.04	C ₄₈ H ₅₀ N ₃ O ₁₃ S ₂ * Na	805-01
Amino-derivative	887.07	C ₄₆ H ₅₄ N ₄ O ₁₀ S ₂	805-02
Maleimide	989.10	C ₅₀ H ₅₃ N ₄ O ₁₂ S ₂ * Na	805-03
			



DY-820

Absorption/emission max.:	823 nm / 843 nm (in Ethanol)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	908.07	C ₄₇ H ₅₄ N ₃ O ₁₀ S ₂ * Na	820-00
NHS-ester	1005.14	C ₅₁ H ₅₇ N ₄ O ₁₂ S ₂ * Na	820-01
Amino-derivative	928.17	C ₄₉ H ₆₁ N ₅ O ₉ S ₂	820-02
Maleimide	1030.19	C ₅₃ H ₆₀ N ₅ O ₁₁ S ₂ * Na	820-03
			



DY-831

Absorption/emission max.:	844 nm / 875 nm (in Ethanol)		
Molar absorbance:	220,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	788.96	C ₄₀ H ₄₉ N ₂ O ₉ S ₂ * Na	831-00
NHS-ester	886.04	C ₄₄ H ₅₂ N ₃ O ₁₁ S ₂ * Na	831-01
Amino-derivative	809.06	C ₄₂ H ₅₆ N ₄ O ₈ S ₂	831-02
Maleimide	911.09	C ₄₆ H ₅₅ N ₄ O ₁₀ S ₂ * Na	831-03
			


DY-845

Absorption/emission max.:	847 nm / 876 nm (in Ethanol)		
Molar absorbance:	160,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1062.19	C ₅₁ H ₅₇ N ₃ O ₁₃ S ₃ * 2 Na	845-00
NHS-ester	1159.26	C ₅₅ H ₆₀ N ₄ O ₁₅ S ₃ * 2 Na	845-01
Amino-derivative	1082.29	C ₅₃ H ₆₄ N ₅ O ₁₂ S ₃ * Na	845-02
Maleimide	1184.31	C ₅₇ H ₆₃ N ₅ O ₁₄ S ₃ * 2 Na	845-03
			

DY-865



Absorption/emission max.:	863 nm / 896 nm (in Ethanol)		
Molar absorbance:	190,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1084.19	C ₅₃ H ₅₅ N ₃ O ₁₃ S ₃ * 2 Na	865-00
NHS-ester	1181.26	C ₅₇ H ₅₈ N ₄ O ₁₅ S ₃ * 2 Na	865-01
Amino-derivative	1104.29	C ₅₅ H ₆₂ N ₅ O ₁₂ S ₃ * Na	865-02
Maleimide	1206.32	C ₅₉ H ₆₁ N ₅ O ₁₄ S ₃ * 2 Na	865-03
			

DYOMICS
Colours for Life





Label | MegaStokes-Dyes for
multi-colour-detection



DY-350XL

Absorption/emission max.:	347 nm / 608 nm (in PBS)		
Molar absorbance:	17,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	634.65	C ₂₅ H ₃₁ N ₄ O ₁₀ S ₂ * Na	350XL-00
NHS-ester	731.74	C ₂₉ H ₃₄ N ₅ O ₁₂ S ₂ * Na	350XL-01
Amino-derivative	654.77	C ₂₇ H ₃₈ N ₆ O ₉ S ₂	350XL-02
Maleimide	756.79	C ₃₁ H ₃₇ N ₆ O ₁₁ S ₂ * Na	350XL-03
			



DY-360XL

Absorption/emission max.:	362 nm / 459 nm (in PBS)		
Molar absorbance:	27,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	632.61	C ₂₇ H ₂₆ N ₂ O ₉ S ₂ * 2 Na	360XL-00
NHS-ester	729.69	C ₃₁ H ₂₉ N ₃ O ₁₁ S ₂ * 2 Na	360XL-01
Amino-derivative	652.71	C ₂₉ H ₃₃ N ₄ O ₈ S ₂ * Na	360XL-02
Maleimide	754.74	C ₃₃ H ₃₂ N ₄ O ₁₀ S ₂ * 2 Na	360XL-03
			



DY-370XL

Absorption/emission max.:	368 nm / 475 nm (in PBS)		
Molar absorbance:	17,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	542.51	C ₂₂ H ₁₉ N ₂ O ₉ S ₂ * Na	370XL-00
NHS-ester	639.59	C ₂₆ H ₂₂ N ₃ O ₁₁ S ₂ * Na	370XL-01
Amino-derivative	562.62	C ₂₄ H ₂₆ N ₄ O ₈ S ₂	370XL-02
Maleimide	664.64	C ₂₈ H ₂₅ N ₄ O ₁₀ S ₂ * Na	370XL-03
			



DY-375XL

Absorption/emission max.:	375 nm / 543 nm (in PBS)		
Molar absorbance:	20,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	447.46	C ₂₀ H ₂₁ N ₃ O ₇ S	375XL-00
NHS-ester	544.53	C ₂₄ H ₂₄ N ₄ O ₉ S	375XL-01
Amino-derivative	526.00	C ₂₂ H ₂₈ N ₅ O ₆ S * Cl	375XL-02
Maleimide	569.58	C ₂₆ H ₂₇ N ₅ O ₈ S	375XL-03
			



DY-376XL

Absorption/emission max.:	387 nm / 549 nm (in PBS)		
Molar absorbance:	25,400 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	576.57	C ₂₃ H ₂₅ N ₂ O ₁₀ S ₂ * Na	376XL-00
NHS-ester	673.64	C ₂₇ H ₂₈ N ₃ O ₁₂ S ₂ * Na	376XL-01
Amino-derivative	596.67	C ₂₅ H ₃₂ N ₄ O ₉ S ₂	376XL-02
Maleimide	698.69	C ₂₉ H ₃₁ N ₄ O ₁₁ S ₂ * Na	376XL-03
			



DY-380XL

Absorption/emission max.:	382 nm / 510 nm (in PBS)		
Molar absorbance:	22,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	448.49	C ₂₁ H ₂₄ N ₂ O ₇ S	380XL-00
NHS-ester	545.56	C ₂₅ H ₂₇ N ₃ O ₉ S	380XL-01
Amino-derivative	527.03	C ₂₃ H ₃₁ N ₄ O ₆ S * Cl	380XL-02
Maleimide	570.61	C ₂₇ H ₃₀ N ₄ O ₈ S	380XL-03
			



DY-395XL

Absorption/emission max.:		397 nm / 572 nm (in PBS)	
Molar absorbance:		22,000 M ⁻¹ cm ⁻¹	
Comments: – soluble in methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	490.53	C ₂₂ H ₂₆ N ₄ O ₇ S	395XL-00
NHS-ester	587.60	C ₂₆ H ₂₉ N ₅ O ₉ S	395XL-01
Amino-derivative	569.07	C ₂₄ H ₃₃ N ₆ O ₆ S * Cl	395XL-02
Maleimide	612.65	C ₂₈ H ₃₂ N ₆ O ₈ S	395XL-03
			



DY-396XL

Absorption/emission max.:		394 nm / 572 nm (in PBS)	
Molar absorbance:		26,600 M ⁻¹ cm ⁻¹	
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	474.53	C ₂₃ H ₂₆ N ₂ O ₇ S	396XL-00
NHS-ester	571.60	C ₂₇ H ₂₉ N ₃ O ₉ S	396XL-01
Amino-derivative	553.07	C ₂₅ H ₃₃ N ₄ O ₆ S * Cl	396XL-02
Maleimide	596.65	C ₂₉ H ₃₂ N ₄ O ₈ S	396XL-03
			



DY-475XL

Absorption/emission max.:	494 nm / 512 nm (in Ethanol)		
Molar absorbance:	90,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	585.62	C ₂₉ H ₃₁ N ₃ O ₁₀ S	475XL-00
NHS-ester	682.69	C ₃₃ H ₃₄ N ₂ O ₁₂ S	475XL-01
Amino-derivative	627.70	C ₃₁ H ₃₇ N ₃ O ₉ S	475XL-02
Maleimide	707.75	C ₃₅ H ₃₇ N ₃ O ₁₁ S	475XL-03
			



DY-480XL

Absorption/emission max.:	504 nm / 631 nm (in Ethanol)		
Molar absorbance:	50,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, DMF, DMSO – bright solid state emission – suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	514.60	C ₂₆ H ₃₀ N ₂ O ₇ S	480XL-00
NHS-ester	611.68	C ₃₀ H ₃₃ N ₃ O ₉ S	480XL-01
Amino-derivative	593.14	C ₂₈ H ₃₇ N ₄ O ₆ S * Cl	480XL-02
Maleimide	636.73	C ₃₂ H ₃₆ N ₄ O ₈ S	480XL-03
			



DY-481XL

Absorption/emission max.:	521 nm / 649 nm (in Ethanol)		
Molar absorbance:	50,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	630.67	C ₂₇ H ₃₁ N ₂ O ₁₀ S ₂ * Na	481XL-00
NHS-ester	727.75	C ₃₁ H ₃₄ N ₃ O ₁₂ S ₂ * Na	481XL-01
Amino-derivative	650.78	C ₂₉ H ₃₈ N ₄ O ₉ S ₂	481XL-02
Maleimide	752.80	C ₃₃ H ₃₇ N ₄ O ₁₁ S ₂ * Na	481XL-03
			



DY-485XL

Absorption/emission max.:	488 nm / 559 nm (in Ethanol)		
Molar absorbance:	50,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in methanol, DMF, DMSO			
– bright solid state emission			
– suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	502.59	C ₂₅ H ₃₀ N ₂ O ₇ S	485XL-00
NHS-ester	599.67	C ₂₉ H ₃₃ N ₃ O ₉ S	485XL-01
Amino-derivative	581.13	C ₂₇ H ₃₇ N ₄ O ₆ S * Cl	485XL-02
Maleimide	624.72	C ₃₁ H ₃₆ N ₄ O ₈ S	485XL-03
			



DY-494XL

Absorption/emission max.:	497 nm / 639 nm (in Ethanol)		
Molar absorbance:	50,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	515.44	C ₂₆ H ₃₁ N ₂ O ₄ · Br	494XL-00
NHS-ester	612.51	C ₃₀ H ₃₄ N ₃ O ₆ · Br	494XL-01
Amino-derivative	549.53	C ₂₈ H ₃₈ N ₄ O ₂ · 2 Cl	494XL-02
Maleimide	637.56	C ₃₂ H ₃₇ N ₄ O ₅ · Br	494XL-03
			



DY-510XL

Absorption/emission max.:	512 nm / 590 nm (in Ethanol)		
Molar absorbance:	50,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in methanol, ethanol, DMF, DMSO – bright solid state emission – suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	554.67	C ₂₉ H ₃₄ N ₂ O ₇ S	510XL-00
NHS-ester	651.74	C ₃₃ H ₃₇ N ₃ O ₉ S	510XL-01
Amino-derivative	633.21	C ₃₁ H ₄₁ N ₄ O ₆ S * Cl	510XL-02
Maleimide	676.80	C ₃₅ H ₄₀ N ₄ O ₈ S	510XL-03
			



DY-511XL

Absorption/emission max.:	510 nm / 595 nm (in PBS)		
Molar absorbance:	50,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	656.70	C ₂₉ H ₃₃ N ₂ O ₁₀ S ₂ * Na	511XL-00
NHS-ester	753.77	C ₃₃ H ₃₆ N ₃ O ₁₂ S ₂ * Na	511XL-01
Amino-derivative	676.80	C ₃₁ H ₄₀ N ₄ O ₉ S ₂	511XL-02
Maleimide	778.82	C ₃₅ H ₃₉ N ₄ O ₁₁ S ₂ * Na	511XL-03
			

DY-520XL

Absorption/emission max.:	522 nm / 662 nm (in Ethanol)		
Molar absorbance:	50,000 M ⁻¹ cm ⁻¹		
Comments:			
– soluble in methanol, DMF, DMSO			
– bright solid state emission			
– suitable for protein labeling, microarray experiments, FisH microscopy, gel electrophoresis			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	514.60	C ₂₆ H ₃₀ N ₂ O ₇ S	520XL-00
NHS-ester	611.68	C ₃₀ H ₃₃ N ₃ O ₉ S	520XL-01
Amino-derivative	593.14	C ₂₈ H ₃₇ N ₄ O ₆ S * Cl	520XL-02
Maleimide	636.73	C ₃₂ H ₃₆ N ₄ O ₈ S	520XL-03
			

DY-521XL


Absorption/emission max.:	526 nm / 666 nm (in Ethanol)		
Molar absorbance:	50,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	630.67	C ₂₇ H ₃₁ N ₂ O ₁₀ S ₂ * Na	521XL-00
NHS-ester	727.75	C ₃₁ H ₃₄ N ₃ O ₁₂ S ₂ * Na	521XL-01
Amino-derivative	650.78	C ₂₉ H ₃₈ N ₄ O ₉ S ₂	521XL-02
Maleimide	752.80	C ₃₃ H ₃₇ N ₄ O ₁₁ S ₂ * Na	521XL-03
			

DYOMICS
Colours for Life




Label | Quencher


DYQ-1

Absorption max.:	543 nm (in PBS)		
Molar absorbance:	44,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	788.75	C ₃₂ H ₃₄ N ₆ O ₁₁ S ₂ * 2 Na	DYQ1-00
NHS-ester	885.83	C ₃₆ H ₃₇ N ₇ O ₁₃ S ₂ * 2 Na	DYQ1-01
Amino-derivative	808.85	C ₃₄ H ₄₁ N ₈ O ₁₀ S ₂ * Na	DYQ1-02
Maleimide	910.88	C ₃₈ H ₄₀ N ₈ O ₁₂ S ₂ * 2 Na	DYQ1-03
Structure on request			


DYQ-2

Absorption max.:	641 nm (in PBS)		
Molar absorbance:	105,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	745.79	C ₃₆ H ₃₆ NO ₁₁ S ₂ * Na	DYQ2-00
NHS-ester	842.86	C ₄₀ H ₃₉ N ₂ O ₁₃ S ₂ * Na	DYQ2-01
Amino-derivative	765.89	C ₃₈ H ₄₃ N ₃ O ₁₀ S ₂	DYQ2-02
Maleimide	867.91	C ₄₂ H ₄₂ N ₃ O ₁₂ S ₂ * Na	DYQ2-03
Structure on request			


DYQ-3

Absorption max.:	683 nm (in PBS)		
Molar absorbance:	80,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	802.84	C ₃₅ H ₃₆ N ₂ O ₁₁ S ₃ * 2 Na	DYQ3-00
NHS-ester	899.91	C ₃₉ H ₃₉ N ₃ O ₁₃ S ₃ * 2 Na	DYQ3-01
Amino-derivative	822.94	C ₃₇ H ₄₃ N ₄ O ₁₀ S ₃ * Na	DYQ3-02
Maleimide	924.97	C ₄₁ H ₄₂ N ₄ O ₁₂ S ₃ * 2 Na	DYQ3-03
Structure on request			


DYQ-4

Absorption max.:	766 nm (in PBS)		
Molar absorbance:	180,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	796.92	C ₄₁ H ₄₅ N ₂ O ₉ S ₂ * Na	DYQ4-00
NHS-ester	893.99	C ₄₅ H ₄₈ N ₃ O ₁₁ S ₂ * Na	DYQ4-01
Amino-derivative	817.02	C ₄₃ H ₅₂ N ₄ O ₈ S ₂	DYQ4-02
Maleimide	919.05	C ₄₇ H ₅₁ N ₄ O ₁₀ S ₂ * Na	DYQ4-03
Structure on request			



DYQ-425

Absorption max.:	429 nm (in PBS)		
Molar absorbance:	24,500 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	784.36	C ₂₃ H ₂₂ Br ₂ N ₄ O ₁₀ S ₂ * 2 Na	425Q-00
NHS-ester	881.43	C ₂₇ H ₂₅ Br ₂ N ₅ O ₁₂ S ₂ * 2 Na	425Q-01
Amino-derivative	804.46	C ₂₅ H ₂₉ Br ₂ N ₆ O ₉ S ₂ * Na	425Q-02
Maleimide	906.48	C ₂₉ H ₂₈ Br ₂ N ₆ O ₁₁ S ₂ * 2 Na	425Q-03
Structure on request			



DYQ-505

Absorption max.:	506 nm (in PBS)		
Molar absorbance:	44,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	713.73	C ₃₁ H ₃₃ N ₅ O ₈ S ₂ * 2 Na	505Q-00
NHS-ester	810.80	C ₃₅ H ₃₆ N ₆ O ₁₀ S ₂ * 2 Na	505Q-01
Amino-derivative	733.83	C ₃₃ H ₄₀ N ₇ O ₇ S ₂ * Na	505Q-02
Maleimide	835.85	C ₃₇ H ₃₉ N ₇ O ₉ S ₂ * 2 Na	505Q-03
Structure on request			


DYQ-660

Absorption max.:	660 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	668.84	C ₃₉ H ₄₄ N ₂ O ₆ S	660Q-00
NHS-ester	765.91	C ₄₃ H ₄₇ N ₃ O ₈ S	660Q-01
Amino-derivative	747.38	C ₄₁ H ₅₁ N ₄ O ₅ S * Cl	660Q-02
Maleimide	790.96	C ₄₅ H ₅₀ N ₄ O ₇ S	660Q-03
			


DYQ-661

Absorption max.:	662 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	770.88	C ₃₉ H ₄₃ N ₂ O ₉ S ₂ * Na	661Q-00
NHS-ester	867.96	C ₄₃ H ₄₆ N ₃ O ₁₁ S ₂ * Na	661Q-01
Amino-derivative	790.99	C ₄₁ H ₅₀ N ₄ O ₈ S ₂	661Q-02
Maleimide	893.01	C ₄₅ H ₄₉ N ₄ O ₁₀ S ₂ * Na	661Q-03
			

DYQ-700


Absorption max.:	691 nm (in PBS)		
Molar absorbance:	58,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	808.89	C ₃₉ H ₄₁ N ₆ O ₈ S ₂ * Na	700Q-00
NHS-ester	905.97	C ₄₃ H ₄₄ N ₇ O ₁₀ S ₂ * Na	700Q-01
Amino-derivative	828.99	C ₄₁ H ₄₈ N ₈ O ₇ S ₂	700Q-02
Maleimide	931.02	C ₄₅ H ₄₇ N ₈ O ₉ S ₂ * Na	700Q-03
Structure on request			

DYOMICS
Colours for Life




PEGylated Fluorescent Labels
reduced nonspecific binding


DY-485XLP4

Absorption/emission max.:	478 nm / 557 nm (in Ethanol)		
Molar absorbance:	50,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	798.89	C ₃₇ H ₅₄ N ₂ O ₁₅ S	485XLP4-00
NHS-ester	895.97	C ₄₁ H ₅₇ N ₃ O ₁₇ S	485XLP4-01
Amino-derivative	877.44	C ₃₉ H ₆₁ N ₄ O ₁₄ S * Cl	485XLP4-02
Maleimide	921.02	C ₄₃ H ₆₀ N ₄ O ₁₆ S	485XLP4-03
Structure on request			


DY-520XLP4

Absorption/emission max.:	514 nm / 665 nm (in Ethanol)		
Molar absorbance:	50,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	838.96	C ₄₀ H ₅₈ N ₂ O ₁₅ S	520XLP4-00
NHS-ester	936.03	C ₄₄ H ₆₁ N ₃ O ₁₇ S	520XLP4-01
Amino-derivative	917.50	C ₄₂ H ₆₅ N ₄ O ₁₄ S * Cl	520XLP4-02
Maleimide	961.08	C ₄₆ H ₆₄ N ₄ O ₁₆ S	520XLP4-03
Structure on request			


DY-547P4

Absorption/emission max.:	557 nm / 571 nm (in PBS)		
Molar absorbance:	150,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1005.17	C ₄₇ H ₆₉ N ₂ O ₁₆ S ₂ * Na	547P4-00
NHS-ester	1102.25	C ₅₁ H ₇₂ N ₃ O ₁₈ S ₂ * Na	547P4-01
Amino-derivative	1025.28	C ₄₉ H ₇₆ N ₄ O ₁₅ S ₂	547P4-02
Maleimide	1127.30	C ₅₃ H ₇₅ N ₄ O ₁₇ S ₂ * Na	547P4-03
Structure on request			


DY-631P4

Absorption/emission max.:	637 nm / 659 nm (in Ethanol)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1061.24	C ₅₀ H ₇₃ N ₂ O ₁₇ S ₂ * Na	631P4-00
NHS-ester	1158.31	C ₅₄ H ₇₆ N ₃ O ₁₉ S ₂ * Na	631P4-01
Amino-derivative	1081.34	C ₅₂ H ₈₀ N ₄ O ₁₆ S ₂	631P4-02
Maleimide	1183.36	C ₅₆ H ₇₉ N ₄ O ₁₈ S ₂ * Na	631P4-03
Structure on request			


DY-647P4

Absorption/emission max.:	656 nm / 674 nm (in PBS)		
Molar absorbance:	250,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1327.53	C ₆₁ H ₉₅ N ₂ O ₂₄ S ₂ * Na	647P4-00
NHS-ester	1424.60	C ₆₅ H ₉₈ N ₃ O ₂₆ S ₂ * Na	647P4-01
Amino-derivative	1347.63	C ₆₃ H ₁₀₂ N ₄ O ₂₃ S ₂	647P4-02
Maleimide	1449.65	C ₆₇ H ₁₀₁ N ₄ O ₂₅ S ₂ * Na	647P4-03
Structure on request			


DY-679P4

Absorption/emission max.:	684 nm / 703 nm (in PBS)		
Molar absorbance:	200,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1631.73	C ₆₉ H ₉₇ N ₂ O ₃₀ S ₄ * 3 Na	679P4-00
NHS-ester	1728.81	C ₇₃ H ₁₀₀ N ₃ O ₃₂ S ₄ * 3 Na	679P4-01
Amino-derivative	1651.83	C ₇₁ H ₁₀₄ N ₄ O ₂₉ S ₄ * 2 Na	679P4-02
Maleimide	1753.86	C ₇₅ H ₁₀₃ N ₄ O ₃₁ S ₄ * 3 Na	679P4-03
Structure on request			


DY-681P4

Absorption/emission max.:	690 nm / 709 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1061.24	C ₅₀ H ₇₃ N ₂ O ₁₇ S ₂ * Na	681P4-00
NHS-ester	1158.31	C ₅₄ H ₇₆ N ₃ O ₁₉ S ₂ * Na	681P4-01
Amino-derivative	1081.34	C ₅₂ H ₈₀ N ₄ O ₁₆ S ₂	681P4-02
Maleimide	1183.36	C ₅₆ H ₇₉ N ₄ O ₁₈ S ₂ * Na	681P4-03
Structure on request			


DY-703P4

Absorption/emission max.:	704 nm / 722 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1125.28	C ₅₄ H ₇₃ N ₂ O ₁₈ S ₂ * Na	703P4-00
NHS-ester	1222.35	C ₅₈ H ₇₆ N ₃ O ₂₀ S ₂ * Na	703P4-01
Amino-derivative	1145.38	C ₅₆ H ₈₀ N ₄ O ₁₇ S ₂	703P4-02
Maleimide	1247.40	C ₆₀ H ₇₉ N ₄ O ₁₉ S ₂ * Na	703P4-03
Structure on request			


DY-704P4

Absorption/emission max.:	705 nm / 721 nm (in Ethanol)		
Molar absorbance:	140,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1227.32	C ₅₄ H ₇₂ N ₂ O ₂₁ S ₃ * 2 Na	704P4-00
NHS-ester	1324.40	C ₅₈ H ₇₅ N ₃ O ₂₃ S ₃ * 2 Na	704P4-01
Amino-derivative	1247.43	C ₅₆ H ₇₉ N ₄ O ₂₀ S ₃ * Na	704P4-02
Maleimide	1349.45	C ₆₀ H ₇₈ N ₄ O ₂₂ S ₃ * 2 Na	704P4-03
Structure on request			


DY-731P4

Absorption/emission max.:	737 nm / 754 nm (in Ethanol)		
Molar absorbance:	220,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1087.27	C ₅₂ H ₇₅ N ₂ O ₁₇ S ₂ * Na	731P4-00
NHS-ester	1184.35	C ₅₆ H ₇₈ N ₃ O ₁₉ S ₂ * Na	731P4-01
Amino-derivative	1107.38	C ₅₄ H ₈₂ N ₄ O ₁₆ S ₂	731P4-02
Maleimide	1209.40	C ₅₈ H ₈₁ N ₄ O ₁₈ S ₂ * Na	731P4-03
Structure on request			


DY-747P4

Absorption/emission max.:	756 nm / 776 nm (in PBS)		
Molar absorbance:	240,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1353.56	C ₆₃ H ₉₇ N ₂ O ₂₄ S ₂ * Na	747P4-00
NHS-ester	1450.64	C ₆₇ H ₁₀₀ N ₃ O ₂₆ S ₂ * Na	747P4-01
Amino-derivative	1373.66	C ₆₅ H ₁₀₄ N ₄ O ₂₃ S ₂	747P4-02
Maleimide	1475.69	C ₆₉ H ₁₀₃ N ₄ O ₂₅ S ₂ * Na	747P4-03
Structure on request			


DY-781P4

Absorption/emission max.:	784 nm / 795 nm (in Ethanol)		
Molar absorbance:	170,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1087.28	C ₅₂ H ₇₅ N ₂ O ₁₇ S ₂ * Na	781P4-00
NHS-ester	1184.35	C ₅₆ H ₇₈ N ₃ O ₁₉ S ₂ * Na	781P4-01
Amino-derivative	1107.38	C ₅₄ H ₈₂ N ₄ O ₁₆ S ₂	781P4-02
Maleimide	1209.40	C ₅₈ H ₈₁ N ₄ O ₁₈ S ₂ * Na	781P4-03
Structure on request			

DY-800P4


Absorption/emission max.:	783 nm / 795 nm (in PBS)		
Molar absorbance:	250,000 M ⁻¹ cm ⁻¹		
Comments: – soluble in water, methanol, ethanol, DMF, DMSO			
Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	1587.77	C ₇₂ H ₁₀₄ N ₂ O ₂₈ S ₃ * 2 Na	800P4-00
NHS-ester	1684.84	C ₇₆ H ₁₀₇ N ₃ O ₃₀ S ₃ * 2 Na	800P4-01
Amino-derivative	1607.87	C ₇₄ H ₁₁₁ N ₄ O ₂₇ S ₃ * Na	800P4-02
Maleimide	1709.89	C ₇₈ H ₁₁₀ N ₄ O ₂₉ S ₃ * 2 Na	800P4-03
Structure on request			

DYOMICS
Colours for Life




Special Dyes

FatRed

Absorption/emission max.:	656 nm / 680 nm (in Ethanol)	
Molar absorbance:	95,000 M ⁻¹ cm ⁻¹	
Comments: – soluble in methanol, ethanol, DMF, DMSO		
Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
697.88	C ₄₁ H ₅₁ N ₃ O ₇	FatRed-00
Structure on request		

Stain01

Absorption/emission max.:	513 nm / 580 nm (in Ethanol)	
Molar absorbance:	120,000 M ⁻¹ cm ⁻¹	
Comments: – soluble in water, methanol, ethanol, DMF, DMSO		
Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
751.94	C ₃₇ H ₅₀ N ₃ O ₈ S ₂ * Na	Stain01-00
Structure on request		



DY-344IN

Absorption/emission max.:	387 nm / 549 nm protonated form at pH 3
Molar absorbance:	27,000 M ⁻¹ cm ⁻¹



Comments:

- soluble in water, methanol, ethanol, DMF, DMSO
- fluorescence indicator, pKs 5.1, bright fluorescence in both forms
- deprotonated form with absorption/ emission max. 342 nm/ 455 nm



Available Modification	Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
Carboxylic acid	454.42	C ₂₀ H ₁₉ N ₂ O ₇ S * Na	344IN-00
NHS-ester	551.50	C ₂₄ H ₂₂ N ₃ O ₉ S * Na	344IN-01
Amino-derivative	474.53	C ₂₂ H ₂₆ N ₄ O ₆ S	344IN-02
Maleimide	576.55	C ₂₆ H ₂₅ N ₄ O ₈ S * Na	344IN-03





V02-07027

Absorption/emission max.:	477 nm / 522 nm (in Ethanol)	
Molar absorbance:	70,000 M ⁻¹ cm ⁻¹	
Comments: – soluble in acetone, ethanol, DMF, DMSO – suitable as DNA-stain		
Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
564.50	C ₂₈ H ₃₇ N ₄ OS · BF ₄	V02-07027
<div style="display: flex; justify-content: space-around; align-items: center;"><div style="text-align: center;"></div><div></div></div>		


V02-08078

Absorption/emission max.:	492 nm / 537 nm (in Ethanol)	
Molar absorbance:	90,000 M-1cm-1	
Comments: – soluble in acetone, ethanol, DMF, DMSO – suitable as DNA-stain		
Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
887.88	C ₃₂ H ₄₆ N ₅ OS * 3 C ₂ F ₃ O ₂	V02-08078
<div><div></div><div></div></div>		



V13-01184

Absorption/emission max.:	481 nm / 526 nm (in Ethanol)	
Molar absorbance:	80,000 M ⁻¹ cm ⁻¹	
Comments: – soluble in acetone, ethanol, DMF, DMSO – suitable as DNA-stain		
Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
578.48	C ₂₈ H ₃₅ N ₄ O ₂ S · BF ₄	V13-01184
<div><div></div><div></div></div>		



MitoDy-1

Absorption max.:	485 nm (in Ethanol)	
Emission max.:	558 nm (in Ethanol)	
Molar absorbance:	62,000 M ⁻¹ cm ⁻¹ (in Ethanol)	
Comments: – nonpolar MegaStokes Dye for multi-colour-detection – Hydroxy-substituted as precursor for phosphoramidite – selective staining of mitochondria		
Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
440.24	C ₂₁ H ₂₅ N ₂ O ₃ + BF ₄	MTD-1
Structure on request		
		



DY-660-SL

Absorption max.:	665 nm (in Ethanol)	
Molar absorbance:	20,000 M ⁻¹ cm ⁻¹	
Comments: – soluble in ethanol, acetonitrile, DMF, DMSO		
Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
620.96	C ₄₃ H ₆₀ N ₂ O	660-SL-00
<div><div></div><div></div></div>		



DY-840-S

Absorption max.:	840 nm (in Ethanol)	
Molar absorbance:	65,000 M ⁻¹ cm ⁻¹	
Comments: – soluble in ethanol, acetonitrile, DMF, DMSO		
Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
559.01	C ₃₂ H ₃₁ N ₂ O · ClO ₄	840-S-00
<div style="display: flex; justify-content: space-around; align-items: center;"><div style="text-align: center;"><p>ClO₄⁻</p></div><div style="text-align: center;"></div></div>		

DY-660-X

Absorption/emission max.:	662 nm / 693 nm (in Dichloromethane)	
Molar absorbance:	80,000 M ⁻¹ cm ⁻¹	
Comments: – soluble in dichloromethane, acetonitrile, DMF, DMSO		
Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
637.09	C ₃₃ H ₃₃ N ₂ O ₅ · ClO ₄	660-X
<div><div></div><div></div></div>		

DY-665-X

Absorption/emission max.:	670 nm / 690 nm (in Dichloromethane)	
Molar absorbance:	80,000 M ⁻¹ cm ⁻¹	
Comments: – soluble in dichloromethane, acetonitrile, DMF, DMSO		
Molecular weight (g · mol ⁻¹)	Molecular formula	Productnumber
516.98	C ₂₆ H ₂₉ N ₂ O ₃ · ClO ₄	665-X
<div><div></div><div></div></div>		

Haptens


Dyomics offers a wide range of labeled haptens like biotin, phalloidin, mycotoxins, amino acids or nucleotides. These haptens can easily be coupled to all of our fluorescent labels.

If you have need for a special hapten please ask us. We will do the coupling with our superior dyes.

Fluorescently labeled biotins

The binding pairs streptavidin – avidin / biotin or antibody / antigen are very common in cellbiology and molecular biology in applications like in situ hybridisation, histochemistry or flow cytometry. With such systems it is possible to probe poor detectable molecules.

The building block for preparing biotin – DY-dye – conjugates is the biotin-PEO3-amine. PEO3 is a spacer consisting of ethyleneglycol units. This spacer between the biotin and the dye separates the biotin moiety from its point of attachment, enhances the ability to bind to the binding sites of avidin or streptavidin and improves the hydrophilic character of the biotin conjugate.




Unit size
500 µg (higher amounts on request)

Fluorescently labeled phalloidins

Phalloidin binds very selectively to actin and allows the generation of impressive images in combination with a fluorophor. 300 units (the standard sales unit) are sufficient for staining 300 microscopy slides. Methanol is an appropriate solvent for the labeled phalloidins.

For all our labels phalloidin-conjugates are available. Please add –33 to the order number of the dye of your choice.



Unit size
10 nmol/300 units (higher amounts on request)

Fluorescently labeled proteins

Fluorescently labeled streptavidin

We offer streptavidin conjugates of all our DY-labels and MegaStokes-Dyes. Each streptavidin conjugate will be delivered with a batch specific absorption spectrum. The conjugates will be shipped in PBS, pH 7.5, 100 mM at a concentration of 1 mg/ml.

Unit size	1 mg at 1 mg/ml (higher amounts on request)
-----------	---

Fluorescently labeled antibodies for secondary detection

Secondary antibodies are very common in applications like fluorescence microscopy and flow cytometry. Dyomics offers several fluorescently labeled polyclonal antibodies for secondary detection, among them

- goat anti-rabbit IgG
- goat anti-mouse IgG &
- goat anti-human IgG

(all affinity purified, H&L). The IgGs will be delivered dissolved in PBS, pH 7.5, 100 mM at a concentration of 1 mg/ml. A batch specific absorption spectrum will be provided.

Unit size	1 mg at 1 mg/ml (higher amounts on request)
-----------	---


Others

If you need an antibody or another protein not listed in this catalogue feel free to contact Dyomics. We can assist you in finding a supplier and we can arrange the import, too. Thus, Dyomics can act as a one stop source for your secondary detection reagents.

Fluorescently labeled nucleotides

For all our labels aminoallyl-conjugates are available. Please add –34 to the order-number of the dye of your choice.

Unit size	100 µl
Concentration	1 mM aqueous solution
Product-number	XXX-34



Protocols

Recommendations / instructions for the labeling of proteins with NHS-esters

Preface:

Fluorescently labelled antibodies are versatile tools in bioanalytics and medical diagnostics. They find their use, among others, in technologies like fluorescence microscopy, flow cytometry, protein microarrays and microtiter plate based detection schemes. There are many suppliers for labelled antibodies, and even more for unconjugated proteins. Kits for the labelling of proteins are also available from various sources. One crucial issue, especially in academic institutions, is the price of the kits or conjugated proteins which is prohibitive for some experiments. The following instructions want to describe an economic method of labelling proteins with simple recipes and instrumentation which allow a full control of all ongoing processes.

The recipes:

Sodium bicarbonate buffer for dialysis and labelling (10x / 500 mM)

Dissolve 21 g (250 mmol) of sodium hydrogen carbonate in 400 ml distilled water. Add 1 g of sodium azide (0.2 per cent, works as bacteriozide). Adjust the pH with a concentrated aqueous solution of sodium hydroxide to pH 8.3. Prior to use, dilute the buffer by adding 9 parts water to one part of the concentrated stock solution (v/v).

100 mM PBS-buffer pH 7.4 for elution and storage

Dissolve 2.19 g sodium dihydrogen phosphate ($\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$, 15.9 mmol), 14.97 g disodium hydrogen phosphate ($\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$, 84.1 mmol), 5.8 g sodium chloride (100 mmol) and 0.32 g sodium azide in 1 l distilled water.

Additional materials and instruments required:

- Protective gloves
- 1 pipet 0.5–10 µl and corresponding tips (white)
- 1 pipet 10–100 µl and corresponding tips (yellow)
- Centrifuge for Eppendorf vials
- Ultrasonification bath (recommended, but not required)
- 1 Pasteur pipet
- Dry DMF (e.g. Fluka PO-no. 40228)
- Sephadex G-25 medium (GE Healthcare)
- Column for size exclusion chromatography (e.g. 5 ml or 10 ml graduated pipettes)
- Silanized glass fiber wadding (e.g. Macherey & Nagel PO-no. 718 002)
- 1 glass bar (at least as long as the column)
- 4 beakers (size can vary)
- 1 tripod with clamp
- Shaker for Eppendorf tubes
- 2 Eppendorf vials (1.5 ml)
- UV-spectrometer

Only when dialysis is required:

- QuixSep Micro dialyzer (Orange Scientific n.v./s.a./ Belgium)
- Dialyzing membrane (also available from Orange Scientific)
- Magnetic stirrer
- Stirring bar

Safety remark:

Always wear protective gloves, lab coat and splash goggles for your personal safety when following the instructions below. In addition, make yourself familiar with the MSDSs which come with the products and chemicals.

Antibody preparation:

Commercially available proteins are usually supplied in aqueous buffer or lyophilised from it. It is essential for a successful labelling reaction that there are no primary or secondary amines present in the reaction solution (otherwise the amine containing buffer components would be labelled instead of the target protein). Thus it is highly recommended to remove buffer components like Tris or ammonia by dialysis against the labelling buffer (50 mM sodium hydrogen carbonate, see recipe above). The QuixSep microdialyzer chamber (see <http://www.orangesci.com/dialqs.htm> for detailed instructions) offers therefore a convenient procedure which is also suitable for low volumes ranging from 100 μ l up to 5 ml. Simply fill your protein solution into the dialysis chamber, place a dialysis membrane which was allowed to equilibrate in the labelling buffer for a few minutes on it and fix it by pushing the collar over it. Dialysis can now be done in a 250 ml beaker against the labelling buffer under constant stirring with a magnetic stirrer. Exchange the labelling buffer twice within three hours. Collect the protein solution from the vial by piercing a pipet through the membrane, aspirating the sample and placing it into an Eppendorf vial. Now you have your protein ready for the labelling reaction.

If you are sure that there are no amines present in your protein solution you can also simply add an appropriate volume of the labelling buffer stock solution to your protein to adjust the pH around 8.3.

To enable a batch-to-batch consistency of the degree of labelling, it is necessary to keep the dye to protein ratio constant for each labelling reaction. The following calculation is an example for the use of DY-647-NHS-ester (MW 761.85 g·mol⁻¹, 1 mg) for labelling 1 mg of an IgG (MW 150,000 g·mol⁻¹). Notwithstanding the fact that an optimal degree of labelling depends on the individual application for the conjugate and can require optimisation, a sixfold molar excess of dye for the labelling reaction is a good value to start with.

The molar ratio between 1 mg dye and 1 mg antibody is 197 (150,000 / 761.85). In order to work with volumes lower than 100 μ l, 98.5 μ l DMF (197 / 2) should be given into the plastic vial containing the 1 mg of DY-647-NHS-ester. To assist proper dissolving of the NHS-ester you can vortex the vial. Alternatively, ultrasonification of the solution for a couple of seconds is also a convenient way to dissolve the reactive dye. It can happen that some solvent and dye is contained in the screw cap of the vial. This material can be recovered by shortly centrifuging the vial.

DMF is superior to the more common DMSO as solvent for reactive dyes for several reasons. It is less hygroscopic, has a lower boiling point (which is important if you want to remove the solvent again from unused dye in a speedvac), and often dyes are better soluble in it. These advantages should override fears about the toxicity of this versatile solvent.

To start the labelling reaction, transfer 3 μ l of the NHS-ester solution to the IgG in the reaction buffer. If you can control the concentration of the protein, a concentration of 5 mg/ml is recommended. In our case, the 1 mg IgG should have been dissolved in 200 μ l of 50 mM bicarbonate buffer. Lower concentrations lead to a reduced labelling efficiency due to the fact that the hydrolysis of the NHS-ester competes the labelling reaction, and the less protein is present the more the hydrolysis dominates.

Allow the labelling to be done in the Eppendorf vial placed in a shaker over the course of two hours. After this time, the NHS-ester should be either bound to the protein or have been hydrolyzed which makes the use of a stopping reagent obsolete.

The Sephadex G-25 medium is a classic gel filtration material. Upon addition of elution buffer, it swells by a factor of five in volume. Thus, it is recommended to put at least 6 ml of elution buffer onto 1.2 g Sephadex in a beaker for filling a 5 ml-column. Double the amounts for a 10 ml column. A 5 ml-column is recommended for reaction solution volumes up to 250 μ l. Please use a 10 ml-column for volumes up to 1 ml.

To prepare the size-exclusion chromatography column, fix the column vertically in the clamp of the tripod. Put a small amount of the silanized glass fiber wadding into the column and push it with the help of the glass bar to the small opening at the ground. The wadding allows the buffer to flow through it while it holds back the gel.

Now pipet once elution buffer into the column, immediately followed by the swollen gel in elution buffer which had been allowed before to equilibrate for at least one hour. Fill the column entirely with the gel and allow thereby the material to settle. If the Sephadex in the beaker runs dry due to the aspiration of the buffer, simply add some more of the elution buffer to liquidize the gel again. Keep a 1 cm region at the top of the column free of the gel.

Now transfer your coloured labelling solution slowly onto the column and allow it to sink into the gel. Wash your pipet carefully (or use a new one) to remove all traces of dye or protein which would contaminate your elution buffer. After the labelling solution has sunk fully into the column, start to elute the conjugate by slowly adding the elution buffer drop by drop onto the column. You can add bigger portions once there is a dye free zone on top of the column. When the reaction solution is eluated, you will see a separation between the labelled protein which runs ahead as a quite sharp band while the free dye is slowly smearing behind.

Once the conjugate arrives at the bottom of the column collect it in an Eppendorf vial. While selecting the size of the vial keep in mind that the volume of the protein conjugate solution increases by a factor of ca. 2.5 on the column. Dispose the used Sephadex gel after collecting the conjugate and wash the column. It would be very time-consuming to wash the free dye from the gel, and traces of it could impurify the next conjugate. It is really not worth the time, and the price for 1 g Sephadex is less than €3.

To characterize your protein conjugate solution it is recommended to record an absorption spectrum of it over the full wavelength range covering both the protein absorption at 280 nm as well as the maximum wavelength of the dye together with the shape of the dye's absorption band. Avoid thereby absorbance values exceeding "2" because such values are hard to interpret. Use cuvettes with smaller pathlengths instead or dilute a small fraction of your conjugate solution appropriately with the elution buffer, e.g. by a factor of 20.

A dye to protein (D/P) ratio can be calculated according to a formula derived from the combination of the Lambert-Beer law for 280 nm and the absorption maximum of the dye. A high uncertainty in this formula is the molar absorbance of the dye at its absorption maximum. Usually, the molar absorbance is determined in organic solvents because this parameter is concentration dependent in aqueous solution.

Thus, one can not be sure if the calculated D/P ratio reflects the experimental truth. The absorption spectrum over the full UV-Vis range is a more reliable analytical information and should be used to compare conjugates among each other.

The 95.5 μ l of the DY-647-NHS-ester which remain from the reaction should be stored at -20°C . An exact shelf life of this solution can not be given. The most crucial factor for the stability of the NHS-ester is to protect it from humidity/moisture. Thus, allow the frozen vial to come completely to room temperature prior to opening it again. A simple test for the remaining activity of the dissolved reactive dye could be to couple it to a small amount (0.5 mg) of a cheap protein like bovine serum albumin (or, in times of BSE, streptavidin) according to the instructions given above when using the dye the first time. Later on, the coupling can be repeated, and the absorption spectrum of the conjugate gives an information about how the NHS-ester degraded over the storage time. In correspondence, the amount of NHS-ester can be adjusted for the labelling reaction with the desired protein.

When labelling proteins with colourless haptens like biotin-NHS-ester, a first coupling with a coloured dye and a cheap protein should be performed on the column whereby the volume of the eluent which leaves the column prior to the conjugate should be measured. Assuming that the elution times for proteins are nearly identical (an allowable assumption since it is about size exclusion chromatography), the same dead volume can be disposed prior to collecting the desired colourless conjugate in a second run on the same column with fresh Sephadex gel.

Services

Aside from our standard product portfolio we offer the customizing of fluorescent labels as well as the synthesis of conjugates on demand. Together with reliable partners we can supply fluorescently labeled DNA-oligonucleotides, fluorogenic, peptide-based enzyme substrates and fluorescently labeled small molecules such as β -estradiol or nucleotides. Feel free to contact Dyomics for a quotation. We are sure that our cost/performance ratio will convince you.

Recommended literature

DY-681 and DY-781 labeled primers

Ritz C.M., Wisseman V., *Microsatellite Analyses of Artificial and Spontaneous Dogrose Hybrids Reveal the Hybridogenic Origin of Rosa micrantha by the Contribution of Unreduced Gametes*, *Journal of Heredity*, doi:10.1093/jhered/esq124.

Review

Napp J., Mathejczyk J.E., Alves F., *Optical imaging in vivo with a focus on paediatric disease: technical progress, current preclinical and clinical applications and future perspectives*, *Pediatr Radiol*, (2011) DOI 10.1007/s00247-010-1907-0.

DY-682-labeled secondary antibodies

Koplin A., Preissler S., Ilina Y., Koch M., Scior A., Erhardt M., Deuerling E., *A dual function for chaperones SSB-RAC and the NAC nascent polypeptide-associated complex on ribosomes*, *The Journal of Cell Biology*, Vol. 189(2010)57–68.

DY-682 and DY-800 labeled secondary antibodies

Ferreira D.N., Katayama I.A., Oliveira I.B., Rosa K.T., Furukawa L.N.S., Coelho M.S., Casarini D.E., Heimann J.C., *Salt-Induced Cardiac Hypertrophy and Interstitial Fibrosis Are Due to a Blood Pressure-Independent Mechanism in Wistar Rats*, *The Journal of Nutrition*, 140(2010)1742-51.

DY-601XL, DY-613, V07-04118 and V07-04146 as stains in acid media

Brockmann S., Arnold T., Schweder B., Bernhard G., *Visualizing Acidophilic Microorganisms in Biofilm Communities Using Acid Stable Fluorescence Dyes*, *Journal of Fluorescence*, 20(2010)943-51.

DY-520XL in TIRF microscopy

Spendier K., Carroll-Portillo A., Lidke K.A., Wilson B.S., Timlin J.A., Thomas J.L., *Distribution and Dynamics of Rat Basophilic Leukemia Immunoglobulin E Receptors (Fc3RI) on Planar Ligand-Presenting Surfaces*, *Biophysical Journal*, 99(2010)388–397

DY-549 and DY-649 on protein microarrays

Schröder C., Jacob A., Tonack S., Radon T.P., Sill M., Zucknick M., Rüffer S., Costello E., Neoptolemos J.P., Crnogorac-Jurcevic T., Bauer A., Fellenberg K., Hoheisel J.D., *Dual-color Proteomic Profiling of Complex Samples with a Microarray of 810 Cancer-related Antibodies*, *Mol Cell Proteomics*, 9(2010)1271–1280.

Review

Paganin-Gioanni A., Bellard E., Paquereau L., Ecochard V., Golzio M., Teissié J., *Fluorescence imaging agents in cancerology*, *Radiol Oncol* 2010; 44(3): 142-148.

DY-520XL in multicolour microscopy

Carroll-Portillo A., Spendier K., Pfeiffer J., Griffiths G., Li H., Lidke K.A., Oliver J.M., Lidke D.S., Thomas J.L., Wilson B.S., Timlin J.A., *Formation of a Mast Cell Synapse: FcRI Membrane Dynamics upon Binding Mobile or Immobilized Ligands on Surfaces*, *The Journal of Immunology*, 184(2010)1328-1338.

DY-752 as imaging agent

Pöschinger T., Janunts E., Brünner H., Langenbucher A., *CCD-based projectional imaging of fluorescent probes in tissue-like media: experimental setup and characterization*, *Zeitschrift für Medizinische Physik*, 2010, in press.

DY-481XL and DYQ-661 as FRET-pair

Tan C., Gajovic-Eichmann N., Stöcklein W.F.M., Polzius R., Bier F.F., Direct detection of 9-tetrahydrocannabinol in saliva using a novel homogeneous competitive immunoassay with fluorescence quenching, *Analytica Chimica Acta*, 658(2010)187-192.

DY-636 on labeled particles

Jiang X., Weise1 S., Hafner M., Röcker C., Zhang F., Parak W.J., Nienhaus G.U., Quantitative analysis of the protein corona on FePt nanoparticles formed by transferrin binding, *J. R. Soc. Interface* 6 February 2010 vol. 7 no. Suppl 1 S5-S13.

DY-682 in tumor imaging

Gong H., Kovar J., Little G., Chen H., Olive D.M., In Vivo Imaging of Xenograft Tumors Using an Epidermal Growth Factor Receptor-Specific Affibody Molecule Labeled with a Near-infrared Fluorophore, *Neoplasia*. 2010 February; 12(2): 139–149.

DY-681-labeled primers

Perovic D., Förster J., Devaux P., Hariri D., Guilleroux M., Kanyuka K., Lyons R., Weyen J., Feuerhelm D., Kastir U., Sourdille P., Röder M., Ordon F., Mapping and diagnostic marker development for Soil-borne cereal mosaic virus resistance in bread wheat, *Mol Breeding*, 23(2009)641–653.

DY-681-UTP

Owczarek-Lipska M., Dolf G., Guziewicz K.E., Leeb T., Schelling C., Posthaus H. Braunschweig M.H., Bovine cardiac troponin I gene (TNNT3) as a candidate gene for bovine dilated cardiomyopathy, *Archiv Tierzucht*, 52(2009)113-123.

DY-781-dUTP in PCR application

Braunschweig M.H., Quantification of global DNA methylation with infrared fluorescence in liver and muscle tissues of differentially fed boars, *Luminescence*, 24(2009)213-216.

DY-633-HaloTag

Kosaka N., Ogawa M., Choyke P.L., Karassina N., Corona C., McDougall M., Lynch D., Hoyt C., Levenson R., Los G.V., Kobayashi H., In Vivo Stable Tumor-Specific Painting in Various Colors Using Dehalogenase-Based Protein-Tag Fluorescent Ligands, *Bioconjug Chem.*, 20(2009)1367–1374.

DY-681 in PCR application

Vollmer T., Kleesiek K., Dreier J., Lipopolysaccharide-binding protein (LBP) gene polymorphisms: Rapid genotyping by real-time PCR and association with infective endocarditis, *Clinical Biochemistry* 42 (2009) 1413–1419.

DY-750 and DY-776 in tumor imaging

Makino A., Kizaka-Kondoh S., Yamahara R., Hara I., Kanzaki T., Ozeki E., Hiraoka M., Kimura S., Near-infrared fluorescence tumor imaging using nanocarrier composed of poly(L-lactic acid)-block-poly(sarcosine) amphiphilic polydepsi-peptide, *Biomaterials*, 30(2009)5156-5160.

DY-485XL labeled secondary antibody

Vag T., Schramm T., Kaiser W.A., Hilger I., Proliferating and quiescent human umbilical vein endothelial cells (HUVECs): a potential in vitro model to evaluate contrast agents for molecular imaging of angiogenesis, *Contrast Media & Molecular Imaging*, 4(2009)192-198.

DY-615 in labeled peptides

Pola R., Studenovsky M., Pechar M., Ulbrich K., Hovorka O., Vetvicka D., Rihova B., HPMa-copolymer conjugates targeted to tumor endothelium using synthetic oligopeptides, *Journal of Drug Targeting*, 2009.

MegaStokes dyes in silica particles

Herz E., Burns A., Bonner D., Wiesner U., Large Stokes-Shift Fluorescent Silica Nanoparticles with ENhanced Emission Over Free Dye for Single Excitation Multiplexing, *Macromolecular Rapid Communications*, 30(2009)1907-10.

DY-635 and DY-647-biotin conjugates

Luschinetz F., Dosche C., Kumke M.U., Influence of Streptavidin on the Absorption and Fluorescence Properties of Cyanine Dyes, *Bioconjugate Chemistry*, 20(2009)576-582.

DY-676, DY-681, DY-731, DY-751 and DY-776 for molecular imaging

Pauli J., Vag T., Haag R., Spieles M., Wenzel M., Kaiser W.A., Resch-Genger U., Hilger I., An in vitro characterization study of new near infrared dyes for molecular imaging, *European Journal of Medicinal Chemistry*, 44(2009)3496-3503.

DY-730, DY-731, DY-732, DY-734 and DY-780 in silica nanoparticles

Herz E., Ow H., Bonner D., Burns A., Wiesner U., Dye structure-optical property correlations in near-infrared fluorescent core-shell silica nanoparticles, *Journal of Materials Chemistry*, DOI:10/1039/b902286d.

V03-01136 (now FatRed) as fat stain

Heilig A., Göggerle A., Hinrichs J., Multiphase visualisation of fat containing beta-lactoglobulin-k-carageenan gels by confocal scanning laser microscopy, using a novel dye, V03-01136, for fat staining, *LWT - Food Science and Technology*, 42(2009)646-53.

DY-547 in antibody labeling

Mindt T.L., Jungi V., Wyss S., Friedli A., Pla G., Novak-Hofer I., Grünberg J., Schibli R., Modification of Different IgG1 Antibodies via Glutamine and Lysine using Bacterial and Human Tissue Transglutaminase, *Bioconjugate Chem.*, 19(2008)271-278.

DY-649 labeled antibodies

Almeida R.S., Brunke S., Albrecht A., Thewes S., Laue M., Edwards J.E.Jr., Filler S.G., Hube B., The Hyphal-Associated Adhesin and Invasin Als3 of *Candida albicans* Mediates Iron Acquisition from Host Ferritin, *PLoS Pathogens*, 4(2008)e1000217.

DY-675 for cancer imaging

Tanisaka H., Kizaka-Kondoh S., Makino A., Tanaka S., Hiraoka M., Kimura S., Near-Infrared Fluorescent Labeled Peptosome for Application to Cancer Imaging, *Bioconjugate Chemistry*, 19(2008)109-117.

DY-485XL in STED-microscopy

Schmidt R., Wurm C.A., Jakobs S., Engelhardt J., Egner A., Hell S.W., Spherical nanosized spot unravels the interior of cells, *nature methods*, 5(2008)539-544.

DY-548 and DY-648 on cDNA-microarrays

Maier M., Wutzler S., Bauer M., Trendafilov P., Henrich D., Marzi I., ALTERED GENE EXPRESSION PATTERNS IN DENDRITIC CELLS AFTER SEVERE TRAUMA: IMPLICATIONS FOR SYSTEMIC INFLAMMATION AND ORGAN INJURY, *SHOCK*, 30(2008)344-351.

DY-676 for in vivo imaging

Deissler V., Rüger R., Frank W., Fahr A., Kaiser W.A., Hilger I., Fluorescent Liposomes as Contrast Agents for In Vivo Optical Imaging of Edemas in Mice, *Small*, 4(2008)1240-6.

DY-676 for in vivo imaging

Lisy M.R., Goermer A., Thomas C., Pauli J., Resch-Genger U., Kaiser W.A., Hilger I., *In Vivo Near-infrared Fluorescence Imaging of Carcinoembryonic Antigen-expressing Tumor Cells in Mice*, *Radiology*, 247(2008)779-787.

DY-505, DY-555 and DY-635 in electrophoresis application

Volke D., Hoffmann R., *Quantitative proteomics by fluorescent labeling of cysteine residues using a set of two cyanine-based or three rhodaminebased dyes*, *ELECTROPHORESIS*, 29(2008)4516-4526.

DY-680 and DY-780 in proteomics

Riederer I.M., Herrero R.M., Leuba G., Riederer B.M., *Serial protein labeling with infrared maleimide dyes to identify cysteine modifications*, *Journal of Proteomics*, 71(2008)222-230.

DY-680 and DY-780 in proteomics (among others) / Review

Riederer B.M., *Non-covalent and covalent protein labeling in two-dimensional gel electrophoresis*, *Journal of Proteomics*, 71(2008)231-244.

DY-547-Streptavidin

Mora J.R., Zielinsky T.L., Nelson B.P., Getts R.C., *Protein detection enhanced by 3DNA dendrimer signal amplification*, *BioTechniques*, 44(2008)815-818.

DY-665-X as standard for NIR-quantum yield measurements

Fischer G.M., Ehlers A.P., Zumbusch A., Daltrozzi E., *Near-Infrared Dyes and Fluorophores Based on Diketopyrrolopyrroles*, *Angew. Chem. Intl. Ed.*, 46(2007)3750-3753.

DY-548 and DY-648 on microarrays:

Auer J., Bläss M., Schulze-Koops H., Russwurm S., Nagel T., Kalden J.R., Röllinghoff M., Beuscher H.U., *Expression and regulation of CCL18 in synovial fluid neutrophils of patients with rheumatoid arthritis*, *Arthritis Research & Therapy*, 9(2007)5R94, 12 pages.

DY-547 on microarrays

Kusnezow W., Banzon V., Schröder C., Schaal R., Hoheisel J.D., Rüffer S., Luft P., Duschl A., Syagailo Y.V., *Antibody microarray-based profiling of complex specimens: systematic evaluation of labeling strategies*, *Proteomics*, 7(2007)1786-1799.

DY-633-labeled secondary antibody

Derwinska K., Sauer U., Preininger C., *Reproducibility of hydrogel slides in on-chip immunoassays with respect to scanning mode, spot circularity, and data filtering*, *Analytical Biochemistry*, 370(2007)38-46.

DY-633-labeled secondary antibody

Derwinska K., Gheber L.A., Preininger C., *A comparative analysis of polyurethane hydrogel for immobilization of IgG on chips*, *Analytica Chimica Acta*, 592(2007)132-138.

DY-485XL in 4Pi microscopy:

Lang M.C., Engelhardt J., Hell S.W., *4Pi microscopy with linear fluorescence excitation*, *Optics Letters*, 32(2007)259-61.

DY-485XL-labeled antibody in microscopy:

Staudt, T., Lang M., Medda R., Engelhardt J., Hell S. W., *2,2'-Thiodiethanol: a new water soluble mounting medium for high resolution optical microscopy*, *Microsc. Res. Tech.*, 70(2007) 1-9.

DY-633-biotin in microscopy:

Rinne J., Albarran B., Jylhävä J., Ihalainen T.O., Kankaanpää P., Hytönen V.P., Stayton P.S., Kulomaa M.S., Vihinen-Ranta M., Internalization of novel non-viral vector TAT-streptavidin into human cells, *BMC Biotechnology*, 7(2007)1-14.

DY-676 as contrast agent:

Lisy M.-R., Hartung A., Lang C., Schüler D., Richter W., Reichenbach J.R., Kaiser W.A., Hilger I., Fluorescent Bacterial Magnetic Nanoparticles as Bimodal Contrast Agents, *Investigative Radiology*, 42(2007)235-41.

DY-547-labeled cGMP:

Biskup C., Kusch J., Schulz E., Nache V., Schwede F., Lehmann F., Hagen V., Benndorf K., Relating ligand binding to activation gating in CNGA2 channels, *nature*, 446(2007)440-3.

DY-415, DY-550 and DY-631 for biochips:

Glaserapp C., Mönch W., Krause H., Zappe H., Biochip reader with dynamic holographic excitation and hyperspectral fluorescence detection, *Journal of Biomedical Optics*, 12(2007)014038, 8 pages.

Volkmer T., Becker C., Prodöhl A., Finger C., Schneider D., Assembly of a transmembrane b-Type cytochrome is mainly driven by transmembrane helix interactions, *Biochimica et Biophysica Acta (BBA) - Biomembranes*, 1758(2006)1815-1822.

DY-676 in in-vivo imaging:

Lisy M.-R., Schüler E., Lehmann F., Czerney P., Kaiser W.A., Hilger I., Diagnosis of peritonitis using near-infrared optical imaging of in vivo labeled monocytes-macrophages, *Journal of Biomedical Optics*, 11(2006)064014, 9 pages.

DY-647-labeled streptavidin:

Kusnezow W., Syagailo Y.V., Rffer S., Baudenstiel N., Gauer C., Hoheisel J.D., Wild D., Goychuk I., Optimal Design of Microarray Immunoassays to Compensate for Kinetic Limitations, *Molecular & Cellular Proteomics*, 5.9(2006)1681-96.

DY-547 and DY-648-labeled benzylguanines:

Kepler A., Arrivoli C., Sironi L., Ellenberg J., Fluorophores for live cell imaging of AGT fusion proteins across the visible spectrum, *BioTechniques*, 41(2006)167-75.

DY-547 and DY-631 in confocal microscopy:

Michalakakis S., Reisert J., Geiger H., Wetzel C., Zong X., Bradley J., Spehr M., Huellett S., Gerstner A., Pfeifer A., Hatt H., Yau K.-W., Biel M., Loss of CNGB1 Protein Leads to Olfactory Dysfunction and Subciliary Cyclic Nucleotide-gated Channel Trapping, *J. Biol. Chem.*, 281(2006)35156-66.

DY-521XL in duplex-RT-PCR:

Dreier J., Stoermer M., Maede D., Burckhardt S., Kleesiek K., Enhanced Reverse Transcription-PCR Assay for Detection of Norovirus Genogroup I, *Journal of Clinical Microbiology*, 44(2006)2714-20.

DY-630 labeled peptide to study GPCRs:

Schneider E., Mayer M., Ziemek R., Li L., Hutzler C., Bernhardt G., Buschauer A., A Simple and Powerful Flow Cytometric Method for the Simultaneous Determination of Multiple Parameters at G Protein-Coupled Receptor Subtypes, *ChemBioChem*, 7(2006)1400-9.

DY-636 in a SNP hairpin probe:

Marmé N., Friedrich A., Denapate D., Hakenbeck R., Knemeyer J.-P., Single nucleotide polymorphism analysis using different colored dye dimer probes, *Chemical Physics Letters*, 428(2006)440-5.

DY-520XL in cross correlation FCS:

Kosturko L.D., Maggipinto M.J., Korza G., Lee J.W., Carson J.H., Barbarese E., *Heterogeneous Nuclear Ribonucleoprotein (hnRNP) E1 Binds to hnRNP A2 and Inhibits Translation of A2 Response Element mRNAs*, *Mol. Biol. Cell*, 17(2006)3521-33.

DY-630-biotin in single molecule detection:

Jung S., Angerer B., Loescher F., Niehren S., Winkle J., Seeger S., *Biotin-Functionalized Cellulose-Based Monolayers as Sensitive Interfaces for the Detection of Single Molecules*, *ChemBioChem*, 7(2006)900-3.

DY-647 and DYQ-660 as FRET pair:

Käppel N.D., Dankbar D.M., Gauglitz G., *DNA Quantification in Nanoliter Volumes*, *Microchim Acta*, 154(2006)65-71.

DY-633-neutravidin as control on peptide microarrays:

Andresen H., Grötzinger C., Zarse K., Kreuzer O.J., Ehrentreich-Förster E., Bier F.F., *Functional peptide microarrays for specific and sensitive antibody diagnostics*, *Proteomics*, 6(2006)1376-84.

DY-633-streptavidin

Andresen H., Zarse K., Grötzinger C., Hollidt J.M., Ehrentreich-Förster E., Bier F.F., Kreuzer O.J., *Development of peptide microarrays for epitope mapping of antibodies against the human TSH receptor*, *Journal of Immunological Methods*, 315(2006)11-18.

DY-781 in sequencing application:

Rothenstein D., Haible D., Dasgupta I., Dutt N., Patil B.L., Jeske H., *Biodiversity and recombination of cassava-infecting begomoviruses from southern India*, *Archives of Virology*, 151(2006)55-69.

DY-681-labeled primers in asymmetric PCR:

de Almeida R.A., Heuser T., Blaschke R., Bartram C.R., Janssen J.W.G., *Control of MYEOV Protein Synthesis by Upstream Open Reading Frames*, *J. Biol. Chem*, 281(2006) 695-704.

DY-633-labeled primers for microarray applications:

Sauer U., Preininger C., Hany-Schmatzberger R., *Quick and simple: quality control of microarray data*, *Bioinformatics*, 21(2005)1572-8.

DY-630 in DNA- and protein arrays:

Preininger C., Sauer U., Dayteg J., Pichler R., *Optimizing processing parameters for signal enhancement of oligonucleotide and protein arrays on ARChip Epoxy*, *Bioelectrochemistry*, 67(2005)155 - 162.

DY-633-labeled peptides:

von Nickisch-Rosenegk M., Ehrentreich-Forster E., Strehlo R., Christmann A., Bier F.F., *Chemically synthesized zinc finger molecules as nano-addressable probes for double-stranded DNAs*, *Journal of Nanobiotechnology*, 3(2005)5.

Detection of DY-681-labeled proteins:

Arndt V., Daniel C., Nastainczyk W., Alberti S., Höhfeld J., *BAG-2 Acts as Inhibitor of the Chaperone-associated Ubiquitin Ligase CHIP*, *Molecular Biology of the Cell*, 16(2005)5891-5900.

DY-547 and DY-631 in microscopy:

Michalakis S., Geiger H., Haverkamp S., Hofmann F., Gerstner A., Biel M., *Impaired Opsin Targeting and Cone Photoreceptor Migration in the Retina of Mice Lacking the Cyclic Nucleotide-Gated Channel CNGA3*, *Investigative Ophthalmology & Visual Science*, 46(2005)1516-1524.

DY-547 in microscopy:

Hüttl S., Michalakakis S., Seeliger, Luo D.-G., Acar N., Geiger H., Hudl K., Mader R., Haverkamp S., Moser M., Pfeifer A., Gerstner A., Yau K.-W., Biel M., *Impaired Channel Targeting and Retinal Degeneration in Mice Lacking the Cyclic Nucleotide-Gated Channel Subunit CNGB1*, *The Journal of Neuroscience*, 25(2005)130-138.

DY-681 and DY-781 in sequencing application:

Kaltenpoth M., Strohm E., Gadau J., *Polymorphic microsatellite markers for a solitary digger wasp, the European beewolf (Philanthus triangulum; Hymenoptera, Sphecidae)*, *Molecular Ecology Notes* 4(2004)589-592.

Review:

Peng X., Draney D., *Near-IR fluorescent dyes for biological applications*, *LabPlus International*, April/May 2004.

EVOblue™ 30 in TIRF-microscopy:

Ruckstuhl T., Seeger S., *Attoliter detection volumes by confocal total-internal-reflection fluorescence microscopy*, *Optics Letters*, 29(2004)568-571.

Detection of DY-630-Maleimid-labeled proteins:




Diekmann S., Siegmund G., Roecker A., Klemm D.O., *Regioselective nitrilotriacetic acid-cellulose-nickel-complexes for immobilisation of His6-tag proteins*, *Cellulose*, 10(2003)53-63.

DY-630 on DNA-microarrays:

Preininger C., Sauer U., *Quality control of chip manufacture and chip analysis using epoxy-chips as model*, *Sensors and Actuators B*, 90(2003)98-103.

Optical properties of DY-630 and DY-635:

Czerney P., Lehmann F., Wenzel M., Buschmann V., Dietrich A., Mohr G.J., *Taylor-Made Dyes for Fluorescence Correlation Spectroscopy (FCS)*, *Biological Chemistry*, 382(2001)495-498.



Our facilities in Jena, Germany.

Dyomics GmbH

Otto-Schott-Str. 15, D-07745 Jena, Germany

Tel.: +49-3641-646 864 · Fax: +49-3641-646 868

e-mail: info@dyomics.com · www.dyomics.com

Design and Layout by:

ctw · gesellschaft für kommunikationsdesign mbH, Jena, Germany

e-mail: info@ctw-jena.de · www.ctw-jena.de

All rights reserved.



Dyomics GmbH

Otto-Schott-Str. 15, D-07745 Jena, Germany
Tel.: +49-3641-646 864 · Fax: +49-3641-646 868
e-mail: info@dyomics.com · www.dyomics.com